NICO BrainPath[®] Instructions for Use


Table of Contents

English	2
Deutsch	10
Español	17
Français	24
Italiano	31
Nederlands	38
Português	45

ENGLISH

NICO BrainPath® Instructions for Use

Caution: Federal (US) Law restricts this device to sale by or on the order of a physician.

Please read all information carefully. Failure to properly follow these instructions may lead to unintended surgical consequences. Important: This package insert is designed to provide instructions for use for the NICO BrainPath. It is not a reference to surgical techniques.

Intended Use

To provide for access and allow for visualization of the surgical field during brain and spinal surgery. Indications may include subcortical access to diseases such as the following:

- Primary/Secondary brain tumors
- Vascular abnormalities/malformations
- Intraventricular tumors/cysts

Contraindications

None known

Warnings

- Use of the NICO BrainPath should be performed only by surgeons having adequate training and experience with using a port to access the brain or spine.
- Avoid use and/or take special care in areas of the brain which carry increased risk.
- Do not resterilize or reuse the NICO BrainPath Sheath or sterile accessories. These sterile devices are intended for a single use only. Resterilization may compromise their integrity, which may result in it malfunctioning.
- The sheath and obturator of the NICO BrainPath which come into contact with bodily fluids require special disposal and handling to prevent biological contamination.
- If prion contamination is known or suspected, it is recommended the reusable components (obturator, manipulation tool and navigation stylet) be incinerated.
- The BrainPath has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of BrainPath in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.


Precautions

Devices and accessories manufactured or distributed by companies not authorized by NICO Corporation may not be compatible with the NICO BrainPath. Use of such products may lead to unanticipated results and possible injury to the user or patient.

Special Facilities, Training, or Qualifications for Users

BrainPath is prescription-use-only and intended to be used by a licensed neurosurgeon. There are no special facility, training, or user qualifications.

Symbols (see primary product label for applicable symbols)

	Prescription Use Only		Catalog Number
	Single use only and should not be reused		Serial Number
	Caution		Non-Sterile
	Indicates that the device is not made with natural rubber latex		Indicates that the device has been sterilized via gamma radiation processes
	Indicates the device is non-pyrogenic		Indicates use-by date
	Indicates Lot Number		Indicates the device manufacturer
	Indicates the manufacture date		Do not use if package is damaged
	Consult the IFU		Indicates item is a Medical Device
	Indicates a double sterile barrier		Unique Device Identifier

Device Description

The BrainPath consists of multiple-sized reusable and re-sterilizable obturators with coordinating single patient use disposable sheaths. The obturator and sheath are assembled in the operating room immediately prior to use. After placement, the obturator is removed leaving behind the sheath which provides a 13.5 mm or 11 mm surgical corridor.

The obturators are designed with a channel and set-screw for placement and securement of third-party navigation probes. The channel of the obturator stops before the tapered tip begins thus navigation settings must reflect the obturator tip length as an extension/offset. For navigation probes that cannot be placed in the existing channel, various adapters are available that ensure the tip of the probe is aligned with the head of the obturator, thus navigation settings must reflect the entire obturator length as an extension/offset. The following tables outline the available BrainPath configurations and relevant obturator dimensions for navigation:


Table 1: BrainPath Configurations

Description	Sheath Inner-Diameter (mm)	Sheath Length (mm)	Obturator Tip Length (mm)
13.5 mm x 50 mm (ST)	13.5	50	7.5
13.5 mm x 50 mm		50	15
13.5 mm x 60 mm		60	
13.5 mm x 75 mm		75	
13.5 mm x 95 mm		95	
11 mm x 50 mm	11	50	
11 mm x 60 mm		60	
11 mm x 75 mm		75	

Table 2: Relevant Obturator Dimensions

Obturator	Third-Party Nav. Probe: Within Obturator Channel		Third-Party Nav. Probe: External using NICO Adapter	
	Acceptable Diameter (mm)	Obturator Tip- Length/ Extension/ Offset (mm)	Acceptable Diameter (mm)	Obturator Overall Length/Extension/ Offset (mm)
13.5 mm x 50 mm (ST)	< 4.3	7.5	4.3 - 5.2	73
13.5 mm x 50 mm		15		81
13.5 mm x 60 mm				91
13.5 mm x 75 mm				106
13.5 mm x 95 mm				126
11 mm x 50 mm				92
11 mm x 60 mm				102
11 mm x 75 mm				116

Figure 1: Obturator Measurement Locations

The BrainPath accessories include a “manipulation tool,” which is similar to a dental probe and is used for manipulating the position of the sheath after it has been placed. The accessories also include “Shepherd’s Hooks” for attaching to various commercially available retractors, a “navigation stylet” for use with the 95 mm obturator, “Probe Adapters” for use with third-party navigation probes that cannot be placed in the channel of the obturator, and a sterilization tray for sterilization of the reusable components (i.e., obturators, manipulation tools, and stylet).

Table 3: NICO BrainPath Components Supplied

BrainPath Component	Sterile or Non-Sterile	Single-Use or Reusable
Sheath (all sizes)	Sterile	Single-Use
Obturators (all sizes)	Non-Sterile	Reusable
Manipulation Tool	Non-Sterile	Reusable
Navigation Stylet	Non-Sterile	Reusable
Sterilization Tray	Non-Sterile	Reusable
Shepherd’s Hooks (all types)	Sterile	Single-Use
Probe Adapter (all types)	Sterile	Single-Use

Instructions for Use

1. The surgeon determines and plans the appropriate approach to the area of interest and prepares the patient according to standard operative preparation protocols as required by the facility.
2. Prior to use, inspect packaging of the sheath and/or sterile accessory (e.g., Shepherd's Hook) for the appropriate size as requested by the surgeon and for any damage or gross deformations. If there is any indication that sterility or functionality has been compromised, obtain a new device and contact NICO Customer Service about returning the original.
3. Remove sheath from packaging using sterile technique.
4. Obtain sterile obturator and confirm that size of obturator matches size of sheath as requested by the surgeon.
5. Obtain sterile manipulation tool and/or navigation stylet (if desired).
6. Obtain appropriate sterile Shepherd's Hook (if desired).
7. Engage the Shepherd's Hook with the proximal portion of the sheath (if desired) just below the knurled ring.
8. Wet the shaft of obturator using sterile saline solution.
9. Wet the entire sheath by dipping it in sterile saline solution
10. Properly assemble obturator and sheath. Ensure that sheath is fully seated on obturator such that there is a smooth transition between the distal end of the sheath and obturator.
11. If desired, place the third-party navigation probe within the obturator channel and secure using thumb screw. Alternatively, for navigation probes that cannot be placed in the obturator channel, obtain the appropriate probe adapter, secure, and then place the third-party navigation probe into the adapter. For the 95 mm Obturator, the 95 mm navigation stylet may be used as a post for attachment of a clamping third-party optical navigation component (if desired).
12. After the creation of an appropriate craniotomy and dural incision, insert BrainPath assembly into tissue to desired location. NOTE: Depth bands exist on obturator and sheath and may be used as a reference during placement.
13. Once the desired location has been reached, advance sheath until the most proximal white band on the obturator appears and is immediately adjacent to the knurled ring on the sheath. Once advanced this distance, the distal end of the sheath will be aligned with the tip of the obturator.
14. While holding the sheath in the desired location, remove obturator from the sheath and the surgical field. Once obturator has been removed, sheath position may be manipulated using manipulation tool. Sheath may also be secured further by placing sutures through the holes present in the knurled ring of the sheath.
15. If using a Shepherd's Hook, secure the proximal end of the hook into the applicable retractor arm.
16. If using a Select Lock Shepherd's Hook, secure the rotational position of the sheath (if desired) by tightening the distal tapered nut with the textured surface.

17. Perform procedure using surgeon desired instrumentation through the BrainPath sheath as appropriate.
18. Once the surgeon has completed the desired end result through the BrainPath, remove sheath and dispose per appropriate hospital procedure.

How Supplied

The BrainPath Sheath and accessories are supplied sterile and packaged for single use. Discard into an appropriate container after use. The BrainPath Obturator, Manipulation Tool and Stylet are packaged non-sterile and require onsite processing/sterilization prior to use. Refer to “BrainPath Reprocessing Instructions” (LL-5904) for detailed information concerning onsite reprocessing.

Special Storage or Handling Conditions

There are no special storage or handling conditions for the BrainPath single-use components. For reusable components, refer to the BrainPath Reprocessing instructions.

Device Related Serious Incident Reporting

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

Clinical Benefits of BrainPath

Published literature on procedures that used BrainPath suggests the following clinical benefits:

- Reduction in parenchymal injury

Clinical Risks

All surgeries carry risks; the following table outlines the adverse events mentioned in the literature inherent to neurosurgical procedures for which BrainPath could be used. These risks are known risks associated with these neurosurgical procedures and underlying diseases. None of the studies stated the events or deaths were attributed to the use of the BrainPath. Post-Market Surveillance data also continues to support low complaints for the BrainPath.

Procedural Risks Table

Event	References	Event	Reference
Aphasia	Mansour et al., 2020	Iatrogenic hypervolemia	Bauer et al., 2017
Arm pain	Marenco-Hillebrand et al., 2020	Ischemia	Mansour et al., 2020
Blood loss	Marenco-Hillebrand et al., 2020	Leg weakness	Marenco-Hillebrand et al., 2020
Brain contusion	Marenco-Hillebrand et al., 2020	MCA perforator	Monroy-Sosa et al., 2021

Cardiac demand	Marenco-Hillebrand et al., 2020	Memory and word recall impairment	Eichber et al., 2020a
Cerebrospinal fluid	Marenco-Hillebrand et al., 2020	Paresthesias	Lin et al., 2020
Death	Bauer et al., 2017; Echeverry et al., 2020; Griessenauer et al., 2018; Przybylowski et al., 2015; Rutkowski et al., 2019; Eliyas et al., 2016;	Postoperative stroke with unknown cause	Marenco-Hillebrand et al., 2020
Deep vein thrombosis	Bakhsheshian et al., 2019; Gassie et al., 2019	Rebleeding	Mansour et al., 2020
Early neurological deficit	Eichberg et al., 2019	Reoperation	Rutkowski et al., 2019
Encephalomyelitis	Gassie et al., 2019	Seizure	Eichber et al., 2020a
Haemorrhage	Mansour et al., 2020	Stroke	Eichberg et al., 2018
Hematoma	Chen et al., 2016	Surgical wound dehiscence	Goren et al., 2019
Hemiparesis	Bakhsheshian et al., 2019; Eliyas et al., 2016	Transient short	Marenco-Hillebrand et al., 2020
Hemiplegia	Akbari et al., 2019; Marenco-Hillebrand et al., 2020; Mansour et al., 2020	Upper extremity	Marenco-Hillebrand et al., 2020
Hemiplegia / hemiparesis	Marenco-Hillebrand et al., 2020	Ventriculitis	Rutkowski et al., 2019
Homonymous / hemianopsia	Lin et al., 2020	Worsening of preoperative dysphasia	Marenco-Hillebrand et al., 2020
Hydrocephalus	Marenco-Hillebrand et al., 2020		

Summary of Safety & Clinical Performance

A summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

Basic-UDIs

Obturator	081674402BP1DA
Sheaths	081674402BP2DC
95 mm Stylet	081674402BP3DE
Shepherd's Hooks	081674402BP4DG
Probe Adapter	081674402BP5DJ
Reusable Components	081674402BP6DL

For More Information

For more information on the NICO BrainPath, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at www.niconeuro.com, including electronic versions of the labeling: www.niconeuro.com/labeling

Patent: www.niconeuro.com/patents

© 2021 NICO Corporation. All Rights Reserved. Printed in the USA.

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© 2021 NICO Corporation. Todos os direitos reservados.