NICO BrainPath® Obturator, Manipulation Tool, and Styler Reprocessing Instructions


Table of Contents

English	2
Deutsch	11
Español	20
Français	29
Italiano	38
Nederlands	47
Português	56

ENGLISH

NICO BrainPath®

Obturator, Manipulation Tool, and Stylet Reprocessing Instructions

Purpose: Proper cleaning and sterilization is critical to achieve effective decontamination and sterilization of medical instrumentation. This manual is intended to assist health care personnel in safe handling, effective reprocessing and maintenance of the NICO BrainPath.

Scope: These reprocessing instructions apply to the obturator, manipulation tool, and stylet components of the NICO BrainPath.

The NICO BrainPath reusable components do not have an indefinite functional life. All reusable instruments are subjected to repeated stresses related to use of the device, cleaning, and sterilization processes.

When using the NICO BrainPath, it is essential that the surgeon and operating personnel are fully conversant with the appropriate technique for the device.

WARNING: Do not resterilize or reuse the NICO BrainPath Sheath. The sheath is intended for a single use only. Resterilization may compromise its integrity, which may result in it malfunctioning.

These reprocessing instructions have been validated as being capable of preparing the NICO BrainPath reusable obturator, manipulation tool, and stylet for reuse. It is the responsibility of the reprocessor to ensure that the reprocessing is actually performed using appropriate equipment, materials, and personnel to achieve the desired result. This normally requires verification and routine monitoring of the process. Any deviation by the reprocessor from these instructions should be evaluated for effectiveness and potential adverse consequences.

Where further information is desired, please contact NICO Corporation Customer Service at **317-660-7118 ext.100**.

Terminology	<ul style="list-style-type: none"> • Chemical – Reactionary substances used for processing/reprocessing (can include detergents, enzymatic cleaners, disinfectants, sterilants or rinse aids). • Cleaning – Process of removing contaminants from an item or surface. • Contaminated – Unsterile or unclean state; made impure by contaminants • Decontamination – Use of physical/chemical means to remove and destroy microorganisms from item or surface; rendering safe for handling or disposal. • Disinfection – A process that destroys pathogens and other microorganisms by physical or chemical means. • Processing/reprocessing – Activity including cleaning, disinfection and sterilization necessary to prepare new or used instrument for its intended use. • Sterilization – A validated process used to render product free from viable microorganisms.
Warnings	<ul style="list-style-type: none"> • If prion contamination is known or suspected, it is recommended the device be incinerated • Follow the instructions and warnings issued by the suppliers of any cleaning and disinfection agents and equipment used. • Rinsing processes should utilize lukewarm water (i.e. water ranging from 22°C to 43°C). • Highly alkaline conditions can damage products with aluminum parts. Alkaline detergent solutions (e.g. Getinge Alkaline, Porlystica Alkaline, etc.) may reduce the use-life of the BrainPath reusable components. Automated cleaning/disinfecting methods should not be used • Complex devices, such as those with tubes, hinges, retractable features, mated surfaces, and textured surface finishes require special attention during cleaning. • Avoid exposure to hypochlorite solutions, as these will promote corrosion. • Scratches or dents can result in breakage. • For instruments produced by another manufacturer, reference the manufacturer’s instructions for use. • Care should be taken to remove any debris, tissue, or bone fragments that may collect on the instrument.
Limitations on Reprocessing	<ul style="list-style-type: none"> • Repeated processing has minimal effects on instrument life and function. • End of useful life is generally determined by wear or damage in surgical use. • Carefully inspect instruments between uses to verify proper functioning. • Damaged instruments should be repaired or replaced to prevent potential patient injury such as loss of metal fragments into the surgical site. • Refer to Inspection & Functional Testing section below.

Reprocessing

<p>Care at the Point of Use</p> <p>¹Cleaning was validated using Enzol® (Cidezyme®) at 1 oz/gal (~30mL/3.8L) at 35-40°C</p>	<ul style="list-style-type: none"> • Clean obturator, manipulation tool, and stylet as soon as possible after use. If cleaning must be delayed, immerse components in enzymatic detergent¹ or water to prevent drying or encrustation of surgical soil. Mix detergent solution per the manufacturer's recommendations using tap water. • Avoid prolonged exposure to saline to minimize the chance of corrosion. • Remove excess soil with a disposable wipe.
<p>Containment and Transportation</p> <p>¹Cleaning was validated using Enzol® (Cidezyme®) at 1 oz/gal (~30mL/3.8L) at 35-40°C</p>	<ul style="list-style-type: none"> • Reprocess components as soon as is reasonably possible or 30 minutes after use. • If reprocessing does not occur within 30 minutes after use, soak components within an enzymatic detergent¹. Mix the solution per the manufacturer's recommendations using tap water. Soak components until reprocessing is possible. • Used instruments must be transported in closed or covered containers to prevent unnecessary contamination risk.
<p>Before First Use</p>	<ul style="list-style-type: none"> • Prior to first use of the BrainPath, the obturators, manipulation tools, and stylet must be cleaned and sterilized per the processes outlined below.
<p>Manual Cleaning</p> <p>¹Cleaning was validated using Enzol® (Cidezyme®) at 1 oz/gal (~30mL/3.8L) at 35-40°C</p>	<ul style="list-style-type: none"> • Thoroughly rinse the obturator, manipulation tool, and stylet under running lukewarm tap water. Pay particular attention to the internal areas. Use a syringe to flush lumens and other hard to reach areas. • Prepare an enzymatic cleaning solution in accordance with the manufacturer's instructions.¹ • Soak soiled components for 20 minutes in the enzymatic solution. Ensure they are fully immersed. • When cleaning, fully immerse components in the cleaning solution to avoid aerosol generation. Brush with a soft non-metallic bristle brush to remove all traces of blood and debris. Pay close attention to threads, crevices, seams, and any hard to reach areas. Use a syringe to flush hard to reach areas. • Remove the instruments from the enzymatic solution and rinse under water for a minimum of 3 minutes. (Thoroughly and aggressively flush difficult to reach areas.) A water jet may facilitate the flushing of instruments. • When rinsing, pay particular attention to internal areas and moveable parts. Use a syringe to flush hard to reach areas. • Dry internal areas with filtered compressed air.

Cleaning Inspection	<ul style="list-style-type: none"> • Inspect obturator, manipulation tool, and stylet before sterilization or storage to ensure the complete removal of soil from surfaces, crevices, and moveable parts. • If areas are difficult to inspect visually, check for blood by immersing or flushing the instrument in a 3% hydrogen peroxide solution. If bubbling is observed, blood is present. • Re-clean components if they are subjected to 3% hydrogen peroxide solution, regardless of whether or not there is evidence of soil/blood.
Disinfection	Obturator, manipulation tool, and stylet must be terminally sterilized prior to surgical use. See sterilization instructions.
Inspection and Functional Testing	<ul style="list-style-type: none"> • Visually inspect the obturator, manipulation tool, and stylet to check for damage and wear. End of life for these components is defined by the following: <ul style="list-style-type: none"> - The presence of any cracks - Surface imperfections such as dents or burrs on the cylindrical portion and/or the conical portion of the obturator - Fading of the markings on the obturator such that they are not visible/legible with an unaided eye - Non-functional set screw on the obturator (is loose or cannot be fully advanced and retracted) - Bend anywhere along the length of the cylindrical portion of the obturator • Visually inspect the sterilization tray and its components for defects. End of life for the tray is defined by the following: <ul style="list-style-type: none"> - Failure of locking mechanism and inability to stay closed - Failure of hinges. - Padded insert deformation to extent that components are not secure throughout the sterilization process - Gross deformation of the tray envelope to extent that functionality (i.e., closing, locking, opening) is compromised.
Packaging	<ul style="list-style-type: none"> • Obturators, manipulation tools, and the stylet must be sterilized in the tray provided by NICO. • All components must be arranged in the tray as indicated in Figures 1 and 2 below to ensure steam penetration to all instrument surfaces. • The tray maybe wrapped in standard medical grade steam sterilization wrap (e.g. Kinguard KC600) using the sequential wrapping technique (refer to ANSI/AAMI ST79). NOTE: Do not perform wrapping for Immediate-Use Steam Sterilization. • For Sterrad sterilization, use a Sterrad approved pouch.

Sterilization	<ul style="list-style-type: none"> • Use a validated, properly maintained and calibrated sterilizer. • Effective sterilization can be achieved using the following cycles: <p><u>Steam Autoclaving Parameters used during Validation</u></p> <table border="0"> <tr> <td>Sterilizer Type:</td> <td>Gravity Displacement</td> <td>Prevacuum</td> </tr> <tr> <td>Configuration:</td> <td>wrapped</td> <td>wrapped</td> </tr> <tr> <td>Temperature:</td> <td>121° C</td> <td>132° C</td> </tr> <tr> <td>Exposure Time:</td> <td>30 minutes</td> <td>4 minutes</td> </tr> <tr> <td>Dry Time:</td> <td>30 minutes</td> <td>30 minutes</td> </tr> </table> <p><u>Recommended Steam Autoclaving Parameters per Annex B of TSE Agents: Safe Working and the Prevention of Infection Guidance</u></p> <table border="0"> <tr> <td>Sterilizer Type:</td> <td>Gravity Displacement</td> <td>Prevacuum</td> </tr> <tr> <td>Configuration:</td> <td>wrapped</td> <td>wrapped</td> </tr> <tr> <td>Temperature:</td> <td>134-137° C</td> <td>134-137° C</td> </tr> <tr> <td>Exposure Time:</td> <td>30 minutes</td> <td>18 minutes</td> </tr> <tr> <td>Dry Time:</td> <td>30 minutes</td> <td>30 minutes</td> </tr> </table> <p><u>*Immediate-Use Steam Sterilization (IUSS)</u></p> <table border="0"> <tr> <td>Sterilizer Type:</td> <td>Gravity Displacement</td> <td>Prevacuum</td> </tr> <tr> <td>Configuration:</td> <td>unwrapped</td> <td>unwrapped</td> </tr> <tr> <td>Temperature:</td> <td>132° C</td> <td>132° C</td> </tr> <tr> <td>Exposure Time:</td> <td>3 minutes</td> <td>3 minutes</td> </tr> </table> <p><u>STERIS V-PRO (H2O2)</u></p> <table border="0"> <tr> <td>Sterilizer:</td> <td>1, 1 Plus, maX, or 60</td> </tr> <tr> <td>Cycles:</td> <td>Non-Lumen, Lumen, or Flexible</td> </tr> <tr> <td>Configuration:</td> <td>double-wrapped</td> </tr> </table> <p>*IUSS is intended for emergency situations such as intraoperative contamination. Instruments processed via IUSS must be used immediately after cooling to room temperature as sterility assurance cannot be maintained for unwrapped items.</p>	Sterilizer Type:	Gravity Displacement	Prevacuum	Configuration:	wrapped	wrapped	Temperature:	121° C	132° C	Exposure Time:	30 minutes	4 minutes	Dry Time:	30 minutes	30 minutes	Sterilizer Type:	Gravity Displacement	Prevacuum	Configuration:	wrapped	wrapped	Temperature:	134-137° C	134-137° C	Exposure Time:	30 minutes	18 minutes	Dry Time:	30 minutes	30 minutes	Sterilizer Type:	Gravity Displacement	Prevacuum	Configuration:	unwrapped	unwrapped	Temperature:	132° C	132° C	Exposure Time:	3 minutes	3 minutes	Sterilizer:	1, 1 Plus, maX, or 60	Cycles:	Non-Lumen, Lumen, or Flexible	Configuration:	double-wrapped
Sterilizer Type:	Gravity Displacement	Prevacuum																																															
Configuration:	wrapped	wrapped																																															
Temperature:	121° C	132° C																																															
Exposure Time:	30 minutes	4 minutes																																															
Dry Time:	30 minutes	30 minutes																																															
Sterilizer Type:	Gravity Displacement	Prevacuum																																															
Configuration:	wrapped	wrapped																																															
Temperature:	134-137° C	134-137° C																																															
Exposure Time:	30 minutes	18 minutes																																															
Dry Time:	30 minutes	30 minutes																																															
Sterilizer Type:	Gravity Displacement	Prevacuum																																															
Configuration:	unwrapped	unwrapped																																															
Temperature:	132° C	132° C																																															
Exposure Time:	3 minutes	3 minutes																																															
Sterilizer:	1, 1 Plus, maX, or 60																																																
Cycles:	Non-Lumen, Lumen, or Flexible																																																
Configuration:	double-wrapped																																																
Storage	Store sterile packaged obturator, manipulation tool, and stylet in a manner that provides protection from dust, moisture, insects, vermin, and extremes of temperature and humidity.																																																

Standard Sets

13.5 mm Set [Blue & Gold]	Qty.	11 mm Set [Green]	Qty.
NN-7026 50mm ST-Gold Obt.	1	--	--
NN-7015 50mm Obt.	1	NN-7036 50mm Obt.	1
NN-7016 60mm Obt.	1	NN-7037 60mm Obt.	1
NN-7017 75mm Obt.	1	NN-7038 75mm Obt.	1
NN-9005 Manipulation Tools	2	NN-9005 Manipulation Tools	2
NN-7035 95mm Obt.	1	--	--
NN-9007 95mm Nav Stylet	1	--	--


Figure 1: 13.5 mm BrainPath Reusable Component Kit


Figure 2: 11 mm BrainPath Reusable Component Kit

Special Storage or Handling Conditions

As noted above, store sterile packaged obturator, manipulation tool, and stylet in a manner that provides protection from dust, moisture, insects, vermin, and extremes of temperature and humidity.

Device Related Serious Incident Reporting

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

Clinical Benefits of BrainPath

Published literature on procedures that used BrainPath suggests the following clinical benefits:

- Reduction in parenchymal injury

Clinical Risks

All surgeries carry risks; the following table outlines the adverse events mentioned in the literature inherent to neurosurgical procedures for which BrainPath could be used. These risks are known risks associated with these neurosurgical procedures and underlying diseases. None of the studies stated the events or deaths were attributed to the use of the BrainPath. Post-Market Surveillance data also continues to support low complaints for the BrainPath.

Event	References	Event	Reference
Aphasia	Mansour et al., 2020	Iatrogenic hypervolemia	Bauer et al., 2017
Arm pain	Marenco-Hillebrand et al., 2020	Ischemia	Mansour et al., 2020
Blood loss	Marenco-Hillebrand et al., 2020	Leg weakness	Marenco-Hillebrand et al., 2020
Brain contusion	Marenco-Hillebrand et al., 2020	MCA perforator	Monroy-Sosa et al., 2021
Cardiac demand	Marenco-Hillebrand et al., 2020	Memory and word recall impairment	Eichber et al., 2020a
Cerebrospinal fluid	Marenco-Hillebrand et al., 2020	Paresthesias	Lin et al., 2020
Death	Bauer et al., 2017; Echeverry et al., 2020; Griessenauer et al., 2018; Przybylowski et al., 2015; Rutkowski et al., 2019; Eliyas et al., 2016;	Postoperative stroke with unknown cause	Marenco-Hillebrand et al., 2020
Deep vein thrombosis	Bakhsheshian et al., 2019; Gassie et al., 2019	Rebleeding	Mansour et al., 2020
Early neurological deficit	Eichberg et al., 2019	Reoperation	Rutkowski et al., 2019
Encephalomyelitis	Gassie et al., 2019	Seizure	Eichber et al., 2020a
Haemorrhage	Mansour et al., 2020	Stroke	Eichberg et al., 2018
Hematoma	Chen et al., 2016	Surgical wound dehiscence	Goren et al., 2019

Hemiparesis	Bakhsheshian et al., 2019; Eliyas et al., 2016	Transient short	Marenco-Hillebrand et al., 2020
Hemiplegia	Akbari et al., 2019; Marenco-Hillebrand et al., 2020; Mansour et al., 2020	Upper extremity	Marenco-Hillebrand et al., 2020
Hemiplegia / hemiparesis	Marenco-Hillebrand et al., 2020	Ventriculitis	Rutkowski et al., 2019
Homonymous / hemianopsia	Lin et al., 2020	Worsening of preoperative dysphasia	Marenco-Hillebrand et al., 2020
Hydrocephalus	Marenco-Hillebrand et al., 2020		

Summary of Safety & Clinical Performance

Where applicable, a summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

Basic-UDIs

Obturators	081674402BP1DA
Sheaths	081674402BP2DC
95 mm Stylet	081674402BP3DE
Shepherd's Hooks	081674402BP4DG
Probe Adapter	081674402BP5DJ
Reusable Components	081674402BP6DL

For More Information

For more information on the NICO BrainPath, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at www.niconeuro.com, including electronic versions of the labeling: www.niconeuro.com/labeling

Patent: www.niconeuro.com/patents

© 2022 NICO Corporation. All Rights Reserved. Printed in the USA.

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© 2022 NICO Corporation. Todos os direitos reservados.