**NICO Myriad-LX™ Illumination Pack**  
**Instructions for Use (IFU)**

Caution: Federal (US) Law restricts this device to sale by or on the order of a physician.

Please read all information carefully. Failure to properly follow these instructions may lead to unintended surgical consequences.

**Intended Use/Purpose**

The Myriad-LX Illumination accessories are intended for use with the Myriad-LX light source for delivery of light to the surgical field to enhance visualization of tissue.

- The light fiber should not be operated for more than 60 minutes continuously
- Refer to the relevant Myriad IFU/Manuals for additional Contraindications, Warnings, and Precautions.


**Contraindications**

None

**Device Description**

The NICO Myriad-LX Illumination Packs consists of two Handpiece Sleeves of different lengths and an Illumination Fiber. The Pack is provided sterile and designed to be utilized only with the NICO Myriad handpiece and the NICO Myriad-LX Light Source. The following table shows the available Illumination Packs

Description	Handpiece Sleeve Options (Qty.)	Fiber Optic (Qty.)
11 GA Illumination Pack	11 GA x 10 cm (1) 11 GA x 13 cm (1)	Illumination Fiber (1)
13 GA Illumination Pack	13 GA x 10 cm (1) 13 GA x 13 cm (1)	

Description	Maximum Concentric Outer Diameter Created	Descriptive Cross-Section Image
11 GA Sleeve	0.241"	
13 GA Sleeve	0.213"	

**Instructions for Use**

NOTE: Refer to Myriad Console and Myriad-LX Light Source Operator's Manuals for additional instructions for setup.

- The appropriate protocol should be employed to prepare the patient for surgery.
- Prior to use, inspect packaging for damage or gross deformations. If there is any indication that sterility or functionality has been compromised, obtain a new Pack and contact NICO Customer Service about returning the original.
- Proper sterile technique should be used to remove the desired length Handpiece Sleeve.
  - NOTE: Two Sleeves are provided of different lengths, but only one is used. Each sleeve features a tag indicating the gauge and length.
- Ensure the proper length Sleeve is selected based on the handpiece that will be used and insert the Myriad Handpiece cannula through the delivery sleeve hub into the clear lumen of the extrusion (Figure 1).

Figure 1: Cannula Passing Through Hub/Clear Lumen

- Snap the Sleeve hub onto the handpiece (Figure 2).

Figure 2: Snapping Delivery Sleeve onto Handpiece

- Connect the proximal tubing clip of the Sleeve onto the tubing exiting the rear of the Handpiece (Figure 3).

Figure 3: Sleeve clip attached onto Handpiece Tubing

- Obtain the Illumination Fiber from the tray and pass the Fiber through the Tuohy-Borst valve of the Sleeve. Advance the Fiber until the depth stop contacts the Tuohy-Borst adapter (Figure 4).

Figure 4: Inserting Fiber into Tuohy-Borst

- Lock the Fiber at the desired location by twisting the Tuohy-Borst clockwise. NOTE: The previous steps can be switched (i.e., Fiber inserted into the Sleeve first then the assembly placed on the handpiece).
- Connect and prime the Handpiece in accordance with the Myriad Handpiece IFU.
- Pass the Illumination Fiber proximal connector (Figure 5) into the non-sterile field for connection to the Myriad-LX Light Source.

Figure 5: Illumination Fiber Proximal Connector

- Power on the Myriad-LX Light Source and adjust to the desired intensity.
  - NOTE: Given the distance to the tissue, low levels of intensity (e.g., < 30%) are likely to provide adequate illumination.
- Once the procedure is completed, dispose of the entire Myriad Handpiece, Sleeve, and Fiber assembly per appropriate hospital procedure.

**Operating Conditions**

Temperature: +15°C to +40°C  
Humidity: 30% to 75%

**Storage and Handling Conditions**

Temperature: -29°C to +60°C  
Humidity: 30% to 80%

**How Supplied**

The Packs are provided sterile for single use. These items are sterilized using ethylene oxide. Discard into an appropriate container after use.

**Device Related Serious Incident Reporting**

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

**Major Risks from Risk Management Activities**

- Electric Shock
- Fire/Explosion
- Damage to Healthy Neurological Tissue
- Excessive Bleeding
- Infection
- Transmission of Disease
- Loss of Function of Equipment in Operating Room due to EMI Emissions
- Tissue Irritation
- Pyrogenic Response (Fever)

The published literature and NICO's post-market complaint processes continue to support very low occurrence or absence of these risks. All residual risks associated with the Myriad continue to be acceptable when weighed against the potential clinical benefits.

**Clinical Benefits**

- Targeted tissue resection
- Favorable safety profile

**Summary of Safety & Clinical Performance**

Where applicable, a summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

**Basic-UDI:** 081674402MY2FW

**For More Information**

For more information on the NICO Myriad, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at [www.niconeuro.com](http://www.niconeuro.com). Additional information may be found at [www.niconeuro.com](http://www.niconeuro.com), including electronic versions of the labelling: [www.niconeuro.com/labelling](http://www.niconeuro.com/labelling)

**Patent:** [www.niconeuro.com/patents](http://www.niconeuro.com/patents)

© 2021 NICO Corporation. All Rights Reserved. Printed in the USA.

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© 2021 NICO Corporation. Reservados todos los derechos. Impreso en los Estados Unidos.