
# myriad | myriad NOVUS™

## Console Operator's Manual

# ENGLISH

# Table of Contents

<b>1</b>	<b>OVERVIEW</b> .....	<b>3</b>
<b>2</b>	<b>WARNINGS AND CAUTIONS</b> .....	<b>4</b>
<b>3</b>	<b>MYRIAD CONSOLE INTRODUCTION</b> .....	<b>6</b>
3.1	COMPONENTS (INCLUDED WITH CONSOLE).....	6
3.2	CONSOLE FUNCTIONS AND SYMBOLS .....	6
3.3	CONSOLE CONNECTIONS .....	10
<b>4</b>	<b>USING THE CONSOLE</b> .....	<b>14</b>
4.1	SYSTEM POWER UP .....	14
4.2	PRIMING THE SYSTEM.....	14
4.3	ASPIRATION AND CUTTING.....	15
<b>5</b>	<b>TROUBLESHOOTING</b> .....	<b>16</b>
5.1	SERVICE AND MAINTENANCE.....	17
5.2	INQUIRIES AND SERVICE .....	17
5.3	MAINTENANCE SCHEDULE .....	18
<b>6</b>	<b>SPECIFICATIONS</b> .....	<b>18</b>
6.1	NICO MYRIAD CONSOLE.....	18
6.2	NICO MYRIAD FOOT PEDAL.....	19
6.3	COMPONENTS & ACCESSORIES LIST .....	20
6.4	CLASSIFICATION .....	21
6.5	CONSOLE EQUIPMENT OPERATING ENVIRONMENTAL CONDITIONS .....	21
6.6	ENVIRONMENTAL .....	21
6.7	SYSTEM ELECTROMAGNETIC EMISSIONS AND IMMUNITY DECLARATIONS.....	21
<b>7</b>	<b>GENERAL INFORMATION</b> .....	<b>24</b>

# 1 Overview

This manual is written for the medical personnel who will be responsible for operating the console. It is extremely important that the operator read and thoroughly understand the contents of this manual and follow the instructions contained herein for reliable, safe and efficient operation of the console.

This manual is for the console and foot pedal components of the system. For more details on the various handpieces or accessories, refer to the Manuals or Instructions for Use (IFUs) provided with them.

## **Indications for Use**

The NICO Myriad™ and Myriad NOVUS™ are powered instruments consisting of a console, handpieces, and accessories intended to perform resection and removal of soft tissue and fluids under direct visualization. Types of direct visualization may include laparoscopic, endoscopic, percutaneous, and open. Applications include those when access to the surgical site is limited, such as Neurosurgical/Spinal and ENT/Otolaryngological.

Specific neurosurgical indications may include diseases such as the following:

- Primary/Secondary Brain Tumors
- Vascular Abnormalities/Malformations (e.g., hemangiomas, cavernomas, and hematoma evacuation)
- Intraventricular Tumors/Cysts

## **Contraindications**

- Gynecologic surgery in which the tissue to be resected is known or suspected to contain malignancy.
- Removal of uterine tissue containing suspected fibroids in patients who are:
  - peri- or post-menopausal, or
  - candidates for en bloc tissue removal (e.g., through vagina or via a mini-laparotomy incision)

**WARNING:** Uterine tissue may contain unsuspected cancer. The use of resection devices during fibroid surgery may spread cancer, and decrease the long-term survival of patients. The information should be shared with patients when considering surgery with the use of these devices.

## **Device Description**

The NICO Myriad™ and Myriad NOVUS™ (hereafter referred to as “Myriad”) are minimally invasive surgical systems specifically designed for the removal of tissues under direct, microscopic, or endoscopic visualization. The design incorporates existing technology and is based on combining minimally invasive, high-speed reciprocating cannula with electronically controlled variable suction in a handpiece capable of precise tissue shaving and rapid tissue debulking. Various handpieces and accessories are available, including the NICO Automated Preservation System™ group of accessories for collection and preservation of resected tissue, and the NOVUS accessories that combine to deliver illumination to the surgical field to aid in visualization.

## 2 Warnings and Cautions

### **CAUTION**

To prevent electric shock, do not use the polarized plug with an extension cord, receptacle or other outlet unless the blades can be fully inserted to prevent blade exposure.

### **WARNINGS**

- To prevent fire or shock hazard, do not expose the console to rain or moisture.
- Use only components with this product that are manufactured by NICO Corporation (NICO) and intended for use with the NICO Myriad System. Use of any other systems' components will void all warranties and may result in system damage.
- To avoid the RISK of electric shock, this equipment must only be connected to a SUPPLY MAINS with protective earth.
- Risk of electric shock - do not open console. To reduce the risk of electric shock, do not remove cover. No user serviceable parts inside. Refer servicing to your distributor or to NICO Customer Service (CS@niconeuro.com).
- No modification of this equipment is allowed.

### **Additional Warnings and Precautions**

- Do not position the equipment such that it is difficult to unplug the power cord from the back of the console.
- The Myriad requires special precautions concerning EMC and needs to be installed and put into service according to the EMC information provided in this document.
- Portable and mobile RF communications equipment can affect the Myriad.
- The Myriad is intended for use by healthcare professionals only. The Myriad may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the Myriad or shielding its location.
- Use of accessories and cables other than those specified, with the exception of accessories and cables qualified and sold by NICO, may result in increased emissions or decreased immunity of the equipment and may cause the Myriad to be non-compliant with the requirements of IEC 60601-1-2
- The Myriad should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- If minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- Do not immerse electrosurgical instruments in liquid.
- Take care not to step on wires and/or tubing connecting console to handpiece during procedures as this may result in compromised functionality and/or unexpected behavior of the Myriad handpiece.
- It is important that the proper collection canister be used to prevent fluid breach into the NICO Myriad Console. The canister should not be overfilled – the level of fluid in the canister should be monitored and the canister should be replaced

when it is full. Only use collection canisters as specified by NICO Corporation. See specifications section of this manual for additional details. In the event of fluid breach into the NICO Myriad Console, discontinue use immediately and contact NICO Corporation to arrange servicing.

- This system should be used only by physicians trained in open or endoscopic tissue removal procedures.
- Products manufactured or distributed by companies not authorized by NICO Corporation may not be compatible with the NICO Myriad System. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all single patient use instruments whether used or unused. Do Not Resterilize any NICO device or accessory if its labeling indicates that it is for single use only. Resterilization of any single use disposable device may compromise the integrity and mechanical function of the instrument or accessory which may result in unintended injury or unacceptable clinical result.
- Should any object or liquid fall into the console, unplug it and have it checked by qualified personnel before further operation.
- Allow adequate air circulation to prevent internal heat build-up. Do not place the console on surfaces or near materials that may prevent proper ventilation.
- Do not install the console in a location near heat sources such as radiators or air ducts, and do not place console in direct sunlight. Isolate console from excessive dust, mechanical vibration or shock.
- Periodically wipe off the console with a soft cloth and a common hospital cleaning solution (including cables if necessary). Stains may be removed with a cloth, lightly dampened with a mild detergent solution. Never use strong solvents or abrasive cleansers as these may damage the enclosure. As a safety precaution, unplug the console before wiping it off.
- Do not place the console on any unstable surfaces which may not properly support it. The console may fall resulting in injury and/or product damage. Use only with a cart, stand or table recommended by the manufacturer or sold with the console. A console and cart combination should be moved with care. Quick stops, excessive force and uneven surfaces may cause the console and cart combination to overturn.
- Grounding reliability can only be achieved when this equipment is connected to an equivalent receptacle marked "Hospital Grade."
- This product is for use only by qualified medical personnel trained for its use.
- **DANGER:** Risk of explosion if used in the presence of flammable anesthetics.
- If the unit has been subjected to a sudden temperature change, moisture may form on the metal inside the unit. If such a temperature change has occurred, allow the unit to achieve room temperature prior turning the unit on.
- Do not damage or modify the power cord. Damage to the power cord may cause a fire or shock hazard. When unplugging the power cord, always hold by the plug and remove it carefully.
- When not in use, unplug the console from the electrical outlet.
- When the console is not being used, store it in a location where it will not be damaged.
- Only ship the console in packaging that has been approved by NICO. Alternate packaging will not guarantee the protection of the device during shipping.
- Consult your distributor or NICO Customer Service if you can't correct a problem using this operator's manual.

- The NICO Myriad System is contraindicated for those patients where increased risk or complications may be associated with the removal or biopsy of tissue based on the physician’s judgment. Patients receiving anticoagulant therapy or who may have bleeding disorders may be at increased risk.
- Precaution should be taken for those patients where increased risk or complications may be associated with the removal or biopsy of tissue based on the physician’s judgment. Patients receiving anticoagulant therapy or who may have bleeding disorders may be at increased risk.
- The Myriad is not designed to coagulate tissue. Prior to resection of highly vascularized abnormalities, devascularization should be performed using alternate instruments.

## 3 Myriad Console Introduction

### 3.1 Components (Included with Console)

- Foot Pedal
- Power Cord
- Nitrogen / Instrument Air Line
- Aspiration Line
- Cart
- Myriad-LX Xenon Light Source (NOVUS version only)

*Refer to the Components and Accessories section in this Manual for a list of Myriad compatible components, including consumables. For additional information, please contact your local sales representative.*


### 3.2 Console Functions and Symbols

#### Front Panel


1. System Power  - Toggle to turn power to the system on or off. I=ON, O=OFF
2. Foot Pedal  - Attachment site for the foot pedal cord to the console.
3. Prime Button  - This button enables the user to prime the handpiece, which

includes simultaneous activation of cutting and aspiration (NOTE: maximum aspiration occurs during priming, regardless of the position of the aspiration control dial). Priming will occur as long as the prime button is depressed. The indicator light above the prime button illuminates while the handpiece is being primed. The cutter enable light, cutter active light and aspiration light will also illuminate while the handpiece is being primed.

4. Aspiration  - Attachment site for the aspiration line. The other end attaches to the top of the single-use collection canister. The aspiration line should always remain attached to the console, except when being replaced or the console is being shipped for service.
5. Aspiration Level Control Dial  - The user can increase or decrease the maximum aspiration level achieved at full depression of the foot pedal using this dial.
6. Aspiration Level Indicator - Digital display indicates aspiration level in non-descript units ranging from 0 to 25.
7. Handpiece  - Attachment site for the handpiece cord to the console. Indicator light above attachment site illuminates during cutting.
8. Cutter Enable Button  - This button enables/disables the cutting mode on the handpiece. The indicator light above this button indicates when cutting is enabled. An audible alert also indicates whether or not cutting is enabled. Two sequential tones indicate that cutting is enabled, while a single tone indicates that cutting is disabled.

## Rear Panel


1. Power Cord Receptacle - Attachment site for detachable hospital grade power cord (IEC320C13).
2. Circuit Breaker - See “Circuit Breaker Reset” in Service and Maintenance section of this manual.
3. Rating Chart - Listing of the electrical specifications for the unit.
4. Duty Cycle Rating – When the cutter is enabled, a duty cycle of 10% is applicable.

5. Nitrogen/Instrument Air Supply Input – Attachment site for the line. The other end of the line will attach to the source. The back of the console is fitted with the male DISS connector. The line is fitted with a female DISS connector at each end. If the OR supply has a different connection, adapters are available for purchase.

## Symbols

	<p>Handpiece - Attachment site for the handpiece cord to the console. Indicator light above attachment site illuminates during cutting.</p>
	<p>Cutter Enable Button - This button enables/disables the cutting mode on the handpiece. The indicator light above this button indicates when cutting is enabled. An audible alert also indicates whether or not cutting is enabled. Two sequential tones indicate that cutting is enabled, while a single tone indicates that cutting is disabled.</p>
	<p>System Power - Toggle to turn power to the system on or off. I=ON, O=OFF</p>
	<p>Foot Pedal - Attachment site for the foot pedal cord to the console.</p>
	<p>Prime Button - This button enables the user to prime the handpiece, which includes simultaneous activation of cutting and aspiration (NOTE: maximum aspiration occurs during priming, regardless of the position of the aspiration control dial). Priming will occur as long as the prime button is depressed. The indicator light above the prime button illuminates while the handpiece is being primed. The cutter enable light, cutter active light and aspiration light will also illuminate while the handpiece is being primed.</p>
	<p>Aspiration - Attachment site for the aspiration hose to the console. The other end of the hose attaches to the collection canister. The indicator light above attachment site illuminates during aspiration.</p>

	<p>Aspiration Level Control Dial - The user can increase or decrease the maximum aspiration level achieved at full depression of the foot pedal using this dial.</p>
	<p>Duty cycle rating symbol - when the cutter is enabled, a duty cycle of 10% is applicable.</p>
	<p>Attention Consult Accompanying Documents</p>
	<p>“Must Refer to Instructions for Use” Symbol</p>
	<p>Interference may occur in the vicinity of known RF transmitting devices and equipment marked with the following symbol.</p>
	<p>Nitrogen/Instrument Supply Input - Attachment site for the nitrogen or instrument air line to the console.</p>
	<p>Symbol for the Waste Electrical and Electronic Equipment Directive (WEEE Directive), European Community directive 2002/96/EC on waste electrical and electronic equipment (WEEE)</p>
	<p>Symbol indicating the catalog number for the device.</p>
	<p>Symbol indicating the serial number for the device.</p>
	<p>Symbol indicating the unique device identifier (UDI) for the device.</p>
	<p>Symbol indicating the location of an applied part. For the Myriad System, the applied part is the metal cannula on the handpiece.</p>
	<p>Symbol which specifies contact information for the European Authorized Representative.</p>
	<p>Symbol indicating the manufacturer of the equipment.</p>

	Symbol indicating the date of manufacture of the equipment.
<b>H<sub>z</sub></b>	Frequency Range
	AC Voltage
<b>A</b>	Current Max
	Circuit Breaker

### 3.3 Console Connections


This section provides information for making all connections necessary to set up the console. The following connections will be covered:

- Power Cord
- Nitrogen / Instrument Air Line
- Aspiration Line
- Canister
- Foot Pedal
- Handpiece

#### Power Cord Connection


The power cord connects the console to the building power supply. Insert the hospital grade power cord into the power cord receptacle on the rear panel (Figure 1). The power cord is to be used for mains disconnection.

Note: Ensure that the power cord is in good condition. A damaged power cord poses an electrical shock hazard. When unplugging the unit, always grasp the plug and pull gently. NEVER pull on the cord.

### **Nitrogen / Instrument Air Line Connection**

Supplied with the console are two approximately 18 ft (~5.5 m) lines with female DISS fittings on each end. One has a label and color code for Instrument Air, while the other is for Nitrogen. Choose proper color line for the hospital supply available. Insert one end of the line into the male fitting on the rear panel of the console (Figure 2). Tighten the connection by rotating the connection nut clockwise. The line should be hand tightened. Do not use tools to tighten the fittings. The opposite end of the hose should be connected to the hospital supply. Before attaching the line to the hospital supply, confirm the supply pressure does not exceed the recommended input supply pressure for the console. A DISS to Schrader adapter is available for Schrader connections. If the hospital supply connection is neither DISS nor Schrader, contact the biomedical engineering department, your distributor, or NICO Customer Service.


### **Aspiration Line Connection**

The console will generate its own aspiration when activated by the foot pedal. During the initial console setup, a reusable flexible aspiration line is connected to the front of the console (Figure 3). To connect, firmly press the end of the tube into the fitting. The tube should insert approximately 0.5 inches (13 mm). The aspiration line should always remain attached to the console, except when being replaced or the console is being shipped for service.

The opposite end of the aspiration line includes a flexible suction connector that is attached to the vacuum port on top of the canister (Figure 4). Ensure that the lid and port caps are secure and sealed to avoid vacuum system leaks. Any leaks in the aspiration system will affect the perform

### **Foot Pedal Connection**

To connect the foot pedal to the console, insert the gray foot pedal connector into the gray receptacle on the console until it clicks into position (Figure 5) (NOTE: the foot pedal connector is keyed and will only connect to the console in one orientation). To remove the connector, pull back on outer sleeve with one hand while securing connector with other hand, then withdraw from console receptacle.

### **Handpiece Connection**

To connect the handpiece to the console, insert the handpiece connector into the receptacle on the console until it clicks into position (Figure 6) (NOTE: the handpiece connector is keyed and will only connect to the console in one orientation). To remove the connector, pull back on the outer sleeve with one hand while securing connector in other hand, then withdraw from console and withdraw from console receptacle.

**Figure 6: Handpiece Connection to Console**

Figure 7 illustrates the handpiece tubing connection directly to the canister lid. The handpiece tubing can also connect to the Tissue Preservation System (TPS) when a specimen collector is being used. Refer to the Handpiece IFU for additional details on potential setup.

If it is desired to collect the tissue removed during the procedure, a specimen collector must be used. Section 5.3 contains a list of accessories and spare parts. Additional information on the various Tissue Preservation System (TPS) setups, which contain the specimen collector, can be found in the Handpiece IFU.

**Figure 7: Handpiece Connection to Canister**


# 4 Using the Console

This section provides information for using the console.


## 4.1 System Power Up

After making all necessary connections, turn the system power on by depressing the switch to the “I” position (Figure 8). At this point, the indicator light next to the power switch will illuminate, other indicator lights on the user face will blink, and a single beep will be heard from the console. The aspiration indicator will read zero.

## 4.2 Priming the System

Before using the system, it is required that each handpiece be primed by aspirating fluid into the aperture, through the handpiece and tubing, and into the canister. This can be accomplished by depressing the prime button (Figure 9) while aspirating sterile fluid from a small bowl or from the reservoir provided in the tray in which the handpiece was packaged. This will ensure that the inside of the handpiece and tubing are lubricated and will reduce the likelihood of tissue occlusions.

### 4.3 Aspiration and Cutting

The console offers two modes of operation: (1) Aspiration only, and (2) Aspiration with Cutting.

This system is supplied with a two-stage foot pedal. The foot pedal consists of a large “main” foot pedal and a smaller “side-kick” foot pedal (Figure 10). When depressed, the large foot pedal will initiate either aspiration only or aspiration with cutting (depending on whether or not cutting has been enabled). The side-kick foot pedal enables the user to toggle between these modes of operation by either enabling or disabling cutting. The user may also toggle between these two modes by depressing the cutter enable button (  ) on the console.

When cutting is enabled, depressing the large foot pedal will result in cutting action as well as aspiration. When cutting is disabled, depressing the large foot pedal will result in aspiration only. The level of aspiration delivered is proportional to the depth to which the foot pedal is depressed. Additionally, the aspiration dial on the console may be adjusted to control the maximum amount of aspiration delivered. A duty cycle of 10% is applicable when cutting is enabled as specified below as well as on the back of the console.

Refer to the Instructions for Use (IFU) provided with the handpiece for additional information concerning the handpiece.

To safely terminate operation of this equipment, cease depression of the foot pedal and power off the console.

# 5 Troubleshooting

This section provides guidance if problems are encountered when operating the console. Please consult the following information before contacting NICO to ensure the problem is not a misunderstanding of the operation of the system. If after reading this section the problem still cannot be resolved, please contact NICO for technical support.

Problem	Possible Cause	Action
No system response	Power cord not plugged in	Check cord. Make sure power cord is plugged in properly.
	Circuit breaker tripped	Reset circuit breaker on back of the console.
No aspiration/suction	Gas source not on	Check to make sure wall source or tank is on and tank is full
Handpiece body becomes warm/hot	Excessive tissue and/or blood in handpiece causing friction	Prime/flush handpiece vigorously with saline/ringers
Handpiece cutting slows down	Excessive tissue and/or blood in handpiece causing friction	Prime/flush handpiece vigorously with saline/ringers
Monitoring equipment in OR is affected by Myriad	Monitoring equipment not compliant with internationally recognized standards for electrical interference	Identify manufacturer and name of monitoring equipment, consult NICO Quality Department
Suction achieved via the foot pedal is less than when pressing the prime button.  Handpiece not aspirating/cutting	Aspiration dial is set too low	Turn the aspiration dial up on the console
	Handpiece not plugged in	Check handpiece connection. Make sure handpiece is plugged into console.
	Aspiration line from console not connected to canister	Check console aspiration line. Make sure it is plugged into the canister.
	Aspiration line from handpiece not connected properly	Check the handpiece aspiration line. Make sure it is plugged into the canister or appropriate TPS accessory.
	Canister leak (crack in canister, crack in lid, uncapped port on lid, loose fit between canister and lid)	Check canister and lid for cracks, make sure lid is on tight, make sure all ports are capped
	Foot pedal not plugged in	Check foot pedal. Make sure foot pedal is plugged into console.
	Aspiration tubing clogged	Prime aspiration tubing. Make sure tubing is free from occlusions.
	Incorrect canister is used – volume is too large for effective vacuum delivery	NICO-specified canister <b>MUST</b> be used with the Myriad console
	Insufficient supply of gas to console – low tank	If a tank is being used, confirm that tank has adequate level. Readout on tank gauge must be at least 80 psi when foot pedal or prime button is fully depressed.
	Insufficient supply of gas to console – input pressure	Gas input pressure from wall or tank must be set at 100 psi.
	Dial setting on console is too low	Turn up aspiration dial on console
	Unvented scope during intraventricular case – all channels of scope are occupied	Explain that one channel of scope must remain open for venting so the handpiece will function properly.

## 5.1 Service and Maintenance

The NICO Myriad System is covered by a one-year warranty for defects, materials and workmanship. Additional extended service contracts are available. Contact your distributor or NICO Corporation directly for details.

**NOTE: Repairs and adjustments inside the console are only to be performed by NICO. Customer attempts to repair/adjust the console will void the initial one-year warranty and/or any extended service contracts.**

### Circuit Breaker Reset

The circuit breaker is located on the rear panel of the console below the power cord receptacle.

1. Unplug console from the power source.
2. Disconnect power cord from receptacle on rear panel.
3. Allow 60 seconds for system to reset prior to proceeding to step 4.
4. Press the button on the circuit breaker to reset.
5. Plug the power cord to the console and into the power source.

### Aspiration Line Replacement

To remove the aspiration line for replacement or shipping:

1. With thumb and forefinger, press the ring at the junction of the aspiration line and aspiration line fitting toward the NICO Myriad console.
2. While the ring is depressed, pull the aspiration line away from the console.

To replace the aspiration line:

1. Firmly press the end of the tube into the fitting. The tube should insert approximately 0.5 inches (13 mm).

## 5.2 Inquiries and Service

Direct inquiries to your distributor or to NICO Corporation:

NICO Corporation  
Customer Service Department  
250 East 96th Street, Suite 125  
Indianapolis, IN 46240

Phone: 317-660-7118 ext.100

If your system requires service, contact your distributor or NICO Corporation Customer Service at 317-660-7118 ext.100. Your distributor or NICO Customer Service will provide you with all necessary information for returning and repairing your console.

## 5.3 Maintenance Schedule

Activity	Occurrence	Action
Wipe off console	As needed	Wipe off console with a clean towel moistened with a hospital approved solution
Wipe off foot pedal	As needed	Wipe off foot pedal with a clean towel moistened with a hospital approved solution
Inspect foot pedal cord	Annually	Check for cuts and damage to outside cover and strain reliefs
Inspect power cord	Annually	Check for cuts and damage to outside cover
Inspect Nitrogen/Instrument Air line	Annually	Check for cuts and damage to outside cover
Test leakage current	Annually	Confirm that leakage current is acceptable per the requirements of the applicable hospital/facility
Inspect markings for legibility	Annually	View markings and confirm they are still legible (free from excessive fading, etc.)

**NOTE:** Of the activities listed above, only leakage current testing must be performed by hospital service personnel.


## 6 Specifications

### 6.1 NICO Myriad Console

Width.....	12.5 in. (318 mm)
Height.....	4.10 in. (104 mm)
Length.....	14.75 in. (375 mm)
Weight.....	16.50 Lbs. (7.0 kg)
Maximum Power.....	40 Watts
AC Voltage Range.....	100-240 V
Frequency Range.....	50-60 Hertz
Max Current.....	0.2 Amps
Circuit Breaker.....	2.0 Amps, 250 V
Nitrogen/Instrument Air Input Pressure Range.....	90-100 PSI (620-690 kPa)
Vacuum Generator.....	0-25.5 inHg (0-65 cmHg) [min. @ sea level]
Power Cord Length.....	15 feet (4.5 m) min.
Nitrogen/Instrument Air Line Length.....	18 feet (5.5 m) min.
Handpiece Cable Length.....	17 feet (5.2 m) min.

## 6.2 NICO Myriad Foot Pedal

Width.....	4.75 in. (121 mm)
Height.....	3.75 in. (95 mm)
Length.....	11.875 in. (302 mm)
Weight.....	6.7 lbs. (3.0 kg)
Cord Length .....	17 feet (5.2 m)

## 6.3 Components & Accessories List

Name	Description	
<b>Capital Components</b>		
Myriad Capital Kit	Cart Console Foot Pedal Nitrogen Line or Instrument Air Line Power Cord	
Myriad NOVUS Capital Kit	Cart Console Myriad-LX Xenon Light Source Foot Pedal Nitrogen Line or Instrument Air Line Power Cord	
<b>Replacement Capital Components</b>		
Nitrogen Line	18 ft. (5.5 m)	
Instrument Air Line	18 ft. (5.5 m)	
Power Cord (US)	15 ft. (5.5 m)	
Power Cord (UK)	15 ft. (5.5 m)	
Aspiration Line	20 in. (508 mm)	
<b>Required Consumable Components</b>		
Canister	800 cc (mL) – <b>NOTE: Use ONLY the canister supplied by NICO.</b>	
Myriad Handpieces	Various Sizes	
<b>Consumable Pack Options</b>		
BrainPath Approach Packs	Various Sizes - includes BrainPath Sheath, Shepherd's Hook, Pre-Bent Myriad Handpiece, meniGLIDE, and SecureHold Retractors	
Illumination Packs	Various Sizes - Myriad-LX Fiber Optic Cable and Myriad-LX Sleeve	
BrainPath Approach + Illumination Pack	Various - Combination of BrainPath Approach and Illumination Packs	
Automated Preservation System Packs	Various – Specimen Collector, Specimen Infusion Valve (SIV), Specimen Preserver, and Specimen Collector Identification Tags	
<b>À La Carte and Optional Accessories</b>		
Specimen Collector & Filter Element - Clamshell	Automated Preservation System (APS) accessories	
Replacement Filter Element - Clamshell		
Specimen Collector & Filter Element – Scoop		
Replacement Filter Element - Scoop		
Specimen Infusion Valve (SIV) with 0.5 mm Metering Line		
Specimen Infusion Valve (SIV) with 0.760 mm Metering Line		
Specimen Preserver		
Handpiece Universal Bending Tool		Handpiece Accessory
Replacement PaediScope Handpiece Adapter		Handpiece Accessory
Replacement LOTTA Handpiece Adapter		Handpiece Accessory
Nitrogen Splitter	Non-Device Accessory – 3-way pneumatic	
DISS/Schrader Adapter	Non-Device Accessory – Non-Swivel Schrader x DISS Male	
DISS/Schrader Swivel Adapter	Non-Device Accessory – Swivel Schrader x DISS Male	
Schrader Quick Connect Adapter	Non- Device Accessory – Male DISS hand-tight / Female Schrader Adapter	

\*Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution company, or customer service for the latest product offerings.

## 6.4 Classification

Class I, Type BF Applied Part  
 Noncontinuous Mode of Operation  
 Console: IPX0 (Ordinary)  
 Foot Pedal: IP68

## 6.5 Console Equipment Operating Environmental Conditions

Temperature range within 15°C to 40°C  
 Relative humidity range within 30% to 75%  
 Atmospheric pressure range within 700 to 1060hPa  
 Operating altitude less than or equal to 3,000 meters above sea level

## 6.6 Environmental

Comply with all local codes when disposing of packaging, equipment or any other portion of this product.

## 6.7 System Electromagnetic Emissions and Immunity Declarations

IEC 60601-1-2:2007 Table 1 Requirements		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The equipment is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the equipment or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## IEC 60601-1-2:2007 Table 2 Requirements

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the Equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 UT = 230 Vac	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment requires continued operation during power mains interruptions, it is recommended that the equipment be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## IEC 60601-1-2:2007 Table 4 Requirements

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the equipment including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup> should be less than the compliance level in each frequency range<sup>b</sup></p> <p>Interference may occur in the vicinity of known RF transmitting devices and equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

IEC 60601-1-2:2007 Table 6 Requirements

Recommended separation distances between portable and mobile RF communications equipment and the Myriad System			
The Myriad is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Myriad can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Myriad as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.24
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38

## 7 General Information

### Device Related Serious Incident Reporting

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

### Summary of Safety & Clinical Performance

Where applicable, a summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

**Basic-UDI:** 081674402MYRIADFX

### For More Information

For more information on the NICO Myriad, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at [www.niconeuro.com](http://www.niconeuro.com).

Patent: [www.niconeuro.com/patents](http://www.niconeuro.com/patents)

© 2020 NICO Corporation. All Rights Reserved. Printed in the USA.

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© 2020 NICO Corporation. Reservados todos los derechos. Impreso en los Estados Unidos.