**NICO Myriad™ Handpiece Instructions for Use**

Please read all information carefully. Failure to properly follow these instructions may lead to unintended surgical consequences. Important: This package insert is designed to provide instructions for use for the NICO Myriad handpiece. It is not a reference to surgical techniques. **Intended Use/Purpose:** For the resection of soft tissue abnormalities found in the central nervous system as determined appropriate by the surgeon.

**Contraindications**

None

**Warnings and Precautions**

**WARNING:** The handpiece and other sterile components are intended for single patient use and should not be reused. Reuse of these devices may result in patient infection and/or device malfunction.

- Removal of tissue using the NICO Myriad System should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications and hazards prior to performance of any minimally invasive procedure.
- This device should only be used by physicians trained in open or endoscopic tissue removal procedures.
- The handpiece should not be operated for more than 60 minutes.
- Take care not to step on wires and/or tubing connecting console to handpiece during procedure as this may result in compromised functionality and/or unexpected behavior of the Myriad handpiece.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by NICO Corporation may not be compatible with the NICO Myriad System. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused. Do not resterilize the handpiece or any other sterile components offered by NICO. These are intended for single patient use only. Resterilization may compromise their integrity, which may result in them malfunctioning.
- The Myriad is not designed to coagulate tissue. Prior to resection of highly vascularized abnormalities, devascularization should be performed using alternate instruments.
- Precaution should be taken for those patients where increased risk or complications may be associated with the removal of tissue based on the physician's judgement. Patients receiving anticoagulant therapy or who may have bleeding disorders may be at increased risk.
- A reduction in aspiration delivery may be due to tissue occlusion within the handpiece. If tissue occlusion is suspected, inspect the handpiece and replace if necessary.
- When using the Myriad handpiece through a neuroendoscope, take care to avoid subjecting the handle portion to surgical fluids exiting the neuroendoscope.
- Flush the handpiece with saline (or the desired solution) between uses during procedures.
- The Myriad handpiece is designed to resect soft tissue. As such, care should be taken by the end user to limit resection to targeted tissue. Unintended resection of tissue may result in injury including but not limited to bleeding. Additionally, care should be taken during resection of vascular abnormalities as bleeding may occur.

**Device Description**

The NICO Myriad System consists of a sterile, single patient use handpiece designed to be utilized only with an approved NICO console and foot pedal. Tissue resection with the handpiece occurs through the combination of aspiration and mechanical cutting action which consists of a non-moving (stationary) outer cannula with a reciprocating inner cutting blade. The handpiece incorporates a distal tissue removal aperture. The handpiece thumbwheel is used to manually turn and orient the tissue removal aperture as desired. The system may also be used in an "aspiration only" mode during which no cutting occurs. Various handpiece sizes are available, including handpieces for potential use in the working channel of certain scopes. The following tables outline the available sizes as well as the potential scope compatibility:

**Standard Handpieces:**

| CANNULA DIAMETER | CANNULA LENGTH | DESCRIPTION   |
|------------------|----------------|---------------|
| 15 gauge         | 10 cm          | 1510          |
|                  | 13 cm          | 1513          |
|                  | 16 cm          | 1516          |
| 13 gauge         | 10 cm          | 1310          |
|                  | 13 cm          | 1313          |
|                  | 16 cm          | Pre-Bent 1313 |
| 11 gauge         | 10 cm          | 1110          |
|                  | 13 cm          | 1113          |
|                  | 16 cm          | Pre-Bent 1113 |

**Working-Channel Handpieces:**

| CANNULA DIAMETER | CANNULA LENGTH | DESCRIPTION              |
|------------------|----------------|--------------------------|
| 19 gauge         | 21.5 cm        | Aesculap PaedScope™      |
|                  | 28 cm          | Karl Storz OI HandyPro™  |
|                  | 28 cm          | Karl Storz Little LOTTA™ |
| 17 gauge         | 21.5 cm        | Karl Storz Devo          |
|                  | 25 cm          | Aesculap MINOP® & InVent |
|                  | 25 cm          | Karl Storz GAAB          |
| 15 gauge         | 25 cm          | Karl Storz LOTTA™        |
|                  | 26.5 cm        | Karl Storz LOTTA™        |

**Instructions for Use – General**

**NOTE:** Refer to Console Operator's Manual for additional instructions for setting up the Myriad System. Refer below for the optional handpiece bending tool and specimen collection instructions.

- The appropriate tool should be employed to prepare the patient for surgery.
- Prior to use, inspect packaging for damage or gross deformations. If there is any indication that sterility or functionality has been compromised, obtain a new handpiece and contact NICO Customer Service about returning the original handpiece.
- Remove handpiece from packaging using sterile technique.
  - Handpiece tray includes a reservoir which may be used for priming if desired.
- Pass handpiece connectors into non-sterile field for connection to console.
- Connect the handpiece blue flexible suction connector to the **front** of canister.
  - Use **ONLY** NICO supplied suction canister.
  - NOTE:** The other blue flexible suction connection coming from the front of the console should be connected to the **top** of the canister as noted in the Myriad Console Operator's Manual.
  - If using the Automated Preservation System™ or Myriad Illumination Pack, the setup may be different – refer to the IFU provided for those accessories.
- Connect the handpiece blue rigid connector to the blue connection port on the front of the console. **[A]**
- To prime handpiece, submerge the aperture at the distal tip of handpiece (**DO NOT SUBMERGE THE ENTIRE HANDPIECE**) into the desired sterile priming fluid and enable prime function on console. Prime handpiece until solution is visible in canister.
- Perform tissue resection using cutting and aspiration as desired. Refer to console Operator's Manual for more information concerning cutting and aspiration functions.
- To change position of aperture during use, rotate dial on top of handpiece. Do **NOT** rotate entire handpiece to reposition aperture.
- Once tissue resection is complete, remove handpiece and dispose per appropriate hospital procedure.

**Instructions for Use – Working-Channel Handpiece**

- Follow the general steps listed above, with the addition of:
  - If using a working channel handpiece (**Figure 1**), a depth control mechanism with luer (**Ia**) or unique MINOP attachment (**Ib**) was added for interfacing with endoscopes.
    - For certain handpieces (e.g., LOTTA, Little LOTTA, and PAEDISCOPE), the depth control mechanism with luer (**Ia**) requires a unique adapter (**2**) to connect to those endoscopes. In addition, for handpieces meant to interface with an OI endoscope, a luer adapter (**3**) is provided to control depth depending on the channel to be used.
    - If applicable, **prior to insertion of handpiece into endoscope**, attach adapter component onto endoscope and secure using thumb screw. Insert the handpiece through adapter component into endoscope and secure via Luer lock.
    - Rotate depth-control dial (**4**) to adjust depth of handpiece during procedure, as clinically appropriate.

Figure 1

**Instructions for Use – Handpiece Bending Tool (optional)**

- Not intended for use on the 25cm versions of the handpieces.**
- Prior to use, inspect packaging for damage or gross deformations. If there is any indication that sterility or functionality has been compromised, obtain a new bending tool and contact NICO Customer Service about returning the original bending tool.
  - Remove bending tool from packaging using sterile technique.
  - If applicable, remove stiffening cannula from handpiece (luer fitting).
  - Align aperture as desired by rotating dial.
  - Slide bending tool over cannula of handpiece. For 11ga handpiece, use the side of the bending tool marked "11ga". For 13ga and 15ga handpieces, use the side of the bending tool marked "13ga, 15ga".
  - Confirm that the wider of the two black bands on the cannula extends past the end of the bending tab prior to bending.
  - (Figure 2)** If this tool does not extend past the end of the tab, the handpiece may not function properly after being bent.
 

**NOTE:** The cannula tip for the 15ga, 10cm and 13ga, 10cm handpieces will not extend past the bending tab as illustrated in **Figure 2**. For these handpieces, simply fully insert them (until they cannot be inserted any further) into the bending tool prior to bending.
  - Bend cannula to desired radius by pressing down on the bending tab with one hand while stabilizing the handpiece with the other hand. Do **NOT** pull down on body of handpiece to bend cannula. The cannula may only be bent one time. Do **NOT** attempt to bend the cannula more than once.
 

**NOTE:** Bending tool **MUST** be used for bending cannula – do **NOT** bend by any other means. Doing so may result in malfunction of the handpiece.
  - Once procedure is complete, discard bending tool per appropriate hospital procedure.


Figure 2

**Operating Conditions**

Temperature: +15°C to +40°C  
Humidity: 30% to 75%

**Storage and Handling Conditions**

Temperature: -29°C to +60°C  
Humidity: 30% to 80%

**How Supplied**

The handpiece, bending tool, and replacement working channel adapters are supplied sterile and packaged for single use. These items are sterilized using gamma irradiation. Discard into an appropriate container after use.

**Device Related Serious Incident Reporting**

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

**Major Risks from Risk Management Activities**

- Electric Shock
- Fire/Explosion
- Damage to Healthy Neurological Tissue
- Excessive Bleeding
- Infection
- Transmission of Disease
- Loss of Function of Equipment in Operating Room due to EMI Emissions
- Tissue Irritation
- Pyrogenic Response (Fever)

The published literature and NICO's post-market complaint processes continue to support very low occurrence or absence of these risks. All residual risks associated with the Myriad continue to be acceptable when weighed against the potential clinical benefits.

**Clinical Benefits:**

- Targeted tissue resection
- Favorable safety profile

**Summary of Safety & Clinical Performance**

Where applicable, a summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

**Basic-UDI:** 081674402MY1FU**For More Information**

For more information on the NICO Myriad, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at [www.niconeuro.com](http://www.niconeuro.com). Additional information may be found at [www.niconeuro.com/labeling](http://www.niconeuro.com/labeling) electronic versions of the labelling: [www.niconeuro.com/labeling](http://www.niconeuro.com/labeling)

**Patent:** [www.niconeuro.com/patents](http://www.niconeuro.com/patents)

© 2024 NICO Corporation. All Rights Reserved. Printed in the USA.

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Atlanta Systems, Ltd.
© 2024 NICO Corporation. Reservados todos los derechos. Impreso en los Estados Unidos.
Atlantic Systems, Ltd.