NICO Myriad Illumination Pack
Instructions for Use (IFU)

Caution: Federal (US) Law restricts this device to sale by or on the order of a physician.

Please read all information carefully. Failure to properly follow these instructions may lead to unintended surgical consequences.

Intended Use/Purpose

The Myriad Illumination Pack is intended for use with the Myriad-LX or Myriad SPECTRA light source for delivery of light to the surgical field.

- The Illumination Fiber should not be operated for more than 60 minutes continuously
- Refer to the relevant Myriad IFU/Manuals for additional Contraindications, Warnings, and Precautions.

Contraindications

None

Warnings and Precautions

- WARNING:** Avoid contact between the distal opening of the delivery sleeve and tissue/debris within the surgical field to prevent complete occlusion. Some surgical approaches (e.g., endonasal) may increase this potential. If unintended reduction in light output occurs and/or occlusion of the sleeve is suspected, both the delivery sleeve and the fiber should be replaced.
- Refer to the relevant Myriad IFU/Manuals for additional Contraindications, Warnings, and Precautions.

Device Description

The NICO Myriad Illumination Pack consists of two Handpiece Sleeves of different lengths and an Illumination Fiber. The Pack is provided sterile and designed to be utilized only with the NICO Myriad handpiece and a NICO light source (LX or SPECTRA). The following table shows the available Illumination Packs:

Description	Handpiece Sleeve Options (Qty.)	Fiber Optic (Qty.)
11 GA Illumination Pack	11 GA x 10 cm (1) 11 GA x 13 cm (1)	Illumination Fiber (1)
13 GA Illumination Pack	13 GA x 10 cm (1) 13 GA x 13 cm (1)	

Description	Maximum Concentric Outer Diameter Created	Descriptive Cross-Section Image
11 GA Handpiece Sleeve	0.241"	
13 GA Handpiece Sleeve	0.213"	

Instructions for Use

NOTE: Refer to Myriad Handpiece, Myriad Console, Myriad-LX Light Source or Myriad SPECTRA Light Source Operator's Manuals for additional instructions for setup.

- The appropriate protocol should be employed to prepare the patient for surgery.
- Prior to use, inspect packaging for damage or gross deformations. If there is any indication that sterility or functionality has been compromised, obtain a new Pack and contact NICO Customer Service about returning the original.
- Proper sterile technique should be used to remove the desired length Handpiece Sleeve.
 - NOTE: Two Handpiece Sleeves are provided of different lengths, but only one is used. Each Handpiece Sleeve features a tag indicating the gauge and length.
- Ensure the proper length Handpiece Sleeve is selected based on the Handpiece that is used and insert the Myriad Handpiece cannula through the Handpiece Sleeve hub into the clear lumen of the extrusion (Figure 1).

Figure 1: Cannula Passing Through Hub/Clear Lumen

- Snap the Handpiece Sleeve hub onto the Handpiece (Figure 2).

Figure 2: Snapping Handpiece Sleeve onto Handpiece

- Connect the proximal tubing clip of the Handpiece Sleeve onto the tubing exiting the rear of the Handpiece (Figure 3).

Figure 3: Handpiece Sleeve clip attached onto Handpiece Tubing

- Obtain the Illumination Fiber from the tray and pass it through the Tuohy-Borst valve of the Handpiece Sleeve. Advance the Illumination Fiber until the depth stop contacts the Tuohy-Borst adapter (Figure 4).

Figure 4: Inserting Illumination Fiber into Tuohy-Borst

- Lock the Illumination Fiber at the desired location by twisting the Tuohy-Borst clockwise. NOTE: The previous steps can be switched (i.e., Illumination Fiber inserted into the Handpiece Sleeve first then the assembly placed on the Handpiece).
- Connect and prime the Handpiece in accordance with the Myriad Handpiece IFU.
- Pass the Illumination Fiber proximal connector (Figure 5) into the non-sterile field for connection to the Myriad-LX Light Source or the Myriad SPECTRA Light Source.


Figure 5: Illumination Fiber Proximal Connector

- Power on the Myriad-LX Light Source or the Myriad SPECTRA Light Source and adjust to the desired intensity.
 - NOTE: When using the Myriad SPECTRA Light Source, additional setup is necessary. Refer to the Myriad SPECTRA Light Source Operator's Manual.
 - NOTE: Given the distance to the tissue, low levels of intensity (e.g., < 30%) for white light are likely to provide adequate illumination.
- Once the procedure is completed, dispose of the entire Myriad Handpiece, Handpiece Sleeve, and Illumination Fiber assembly per appropriate hospital procedure.

Operating Conditions

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distributor/ Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

Temperature: +15°C to +40°C

Humidity: 30% to 75%

Storage and Handling Conditions

Temperature: -29°C to +60°C

Humidity: 30% to 80%

How Supplied

The NICO Myriad Illumination Pack is provided sterile for single use. It is sterilized using ethylene oxide. Discard into an appropriate container after use.

Device Related Serious Incident Reporting

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distributor/ Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

Major Risks from Risk Management Activities

- Electric Shock
- Fire/Explosion
- Damage to Healthy Neurological Tissue
- Excessive Bleeding
- Infection
- Transmission of Disease
- Loss of Function of Equipment in Operating Room due to EMI Emissions
- Tissue Irritation
- Pyrogenic Response (Fever)

The published literature and NICO's post-market complaint processes continue to support very low occurrence or absence of these risks. All residual risks associated with the Myriad continue to be acceptable when weighed against the potential clinical benefits.

Clinical Benefits:

- Targeted tissue resection
- Favorable safety profile

Summary of Safety & Clinical Performance

Where applicable, a summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

Basic-UDI: 081674402MY2FW

Symbol

Symbol indicating a Type BF applied part. The Illumination Fiber is the applied part for the Light Source.

For More Information

For more information on the NICO Myriad, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at www.niconero.com. Additional information may be found at www.niconero.com, including electronic versions of the labelling: www.niconero.com/labelling

Patent: www.niconero.com/patents

© 2024 NICO Corporation. All Rights Reserved. Printed in the USA.

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