myriad | myriad
NOVUS™

MYRIAD
SPECTRA™

Console
Operator's Manual



Table of Contents

ENGLISH	1
NEDERLANDS	27
FRANCAIS.....	53
DEUTSCH	79
ITALIANO.....	105
PORTUGUES	131
ESPANOL	157

ENGLISH

Table of Contents

- 1 OVERVIEW.....3
- 2 WARNINGS AND CAUTIONS3
- 3 MYRIAD CONSOLE INTRODUCTION5
 - 3.1 COMPONENTS (INCLUDED WITH CONSOLE)..... 5
 - 3.2 CONSOLE FUNCTIONS AND SYMBOLS 6
 - 3.3 CONSOLE CONNECTIONS 10
- 4 USING THE CONSOLE13
 - 4.1 SYSTEM POWER UP 13
 - 4.2 PRIMING THE SYSTEM 14
 - 4.3 ASPIRATION AND CUTTING..... 14
- 5 TROUBLESHOOTING.....15
 - 5.1 SERVICE AND MAINTENANCE..... 16
 - 5.2 INQUIRIES AND SERVICE 17
 - 5.3 MAINTENANCE SCHEDULE 17
- 6 SPECIFICATIONS18
 - 6.1 NICO MYRIAD AND MYRIAD NOVUS CONSOLES..... 18
 - 6.2 NICO MYRIAD SPECTRA CONSOLE 18
 - 6.3 NICO MYRIAD FOOT PEDAL..... 19
 - 6.4 CLASSIFICATION 19
 - 6.5 CONSOLE EQUIPMENT OPERATING ENVIRONMENTAL CONDITIONS 19
 - 6.6 DISPOSAL..... 19
 - 6.7 COMPONENTS & ACCESSORIES LIST 20
 - 6.8 SYSTEM ELECTROMAGNETIC EMISSIONS AND IMMUNITY DECLARATIONS..... 21
- 7 GENERAL INFORMATION25

1 Overview

This manual is written for the medical personnel who will be responsible for operating the console. It is extremely important that the operator read and thoroughly understand the contents of this manual and follow the instructions contained herein for reliable, safe and efficient operation of the console.

This manual is for the console and foot pedal components of the system. For more details on the various handpieces or accessories, refer to the Manuals or Instructions for Use (IFUs) provided with them.

Intended Use/Purpose

For the resection of soft tissue abnormalities found in the central nervous system as determined appropriate by the surgeon.

Contraindications

None

Special Facilities, Training, or Qualifications for Users

Myriad is prescription-use-only and intended to be used by a licensed surgeon. There are no special facility, training, or user qualifications.

Device Description

The NICO Myriad[®], Myriad NOVUS[®] and Myriad SPECTRA[™] (hereafter referred to as “Myriad”) are minimally invasive surgical systems specifically designed for the removal of tissues under direct, microscopic, or endoscopic visualization. The design incorporates existing technology and is based on combining minimally invasive, high-speed reciprocating cannula with electronically controlled variable suction in a handpiece capable of precise tissue shaving and rapid tissue debulking. Various handpieces and accessories are available, including the NICO Automated Preservation System[™] group of accessories for collection and preservation of resected tissue, and the illumination accessories that combine to deliver illumination to the surgical field to aid in visualization.

2 Warnings and Cautions

CAUTION

To prevent electric shock, do not use the polarized plug with an extension cord, receptacle or other outlet unless the blades can be fully inserted to prevent blade exposure.

WARNINGS

- To prevent fire or shock hazard, do not expose the console to rain or moisture.
- Use only components with this product that are manufactured by NICO Corporation (NICO) and intended for use with the NICO Myriad System. Use of any other systems' components will void all warranties and may result in system damage.
- To avoid the RISK of electric shock, this equipment must only be connected to a SUPPLY MAINS with protective earth.
- Risk of electric shock - do not open console. To reduce the risk of electric shock, do not remove cover. No user serviceable parts inside. Refer servicing to your distributor or to NICO Customer Service (CS@niconeuro.com).
- No modification of this equipment is allowed.

Additional Warnings and Precautions

- Do not position the equipment such that it is difficult to unplug the power cord from the back of the console.
- The Myriad requires special precautions concerning EMC and needs to be installed and put into service according to the EMC information provided in this document.
- Portable and mobile RF communications equipment can affect the Myriad.
- The Myriad is intended for use by healthcare professionals only. The Myriad may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the Myriad or shielding its location.
- Use of accessories and cables other than those specified, with the exception of accessories and cables qualified and sold by NICO, may result in increased emissions or decreased immunity of the equipment and may cause the Myriad to be non-compliant with the requirements of IEC 60601-1-2.
- The Myriad should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- If minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- Do not immerse electrosurgical instruments in liquid.
- Take care not to step on wires and/or tubing connecting console to handpiece during procedures as this may result in compromised functionality and/or unexpected behavior of the Myriad handpiece.
- It is important that the proper collection canister be used to prevent fluid breach into the Console. The canister should not be overfilled – the level of fluid in the canister should be monitored and the canister should be replaced when it is full. Only use collection canisters as specified by NICO Corporation. See specifications section of this manual for additional details. In the event of fluid breach into the Console, discontinue use immediately and contact NICO Corporation to arrange servicing.
- This system should be used only by physicians trained in open or endoscopic tissue removal procedures.
- Products manufactured or distributed by companies not authorized by NICO Corporation may not be compatible with the Myriad. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all single patient use instruments whether used or unused. Do Not Resterilize any NICO device or accessory if its labeling indicates that it is for single use only. Resterilization of any single use disposable device may compromise the integrity and mechanical function of the instrument or accessory which may result in unintended injury or unacceptable clinical result.
- Should any object or liquid fall into the console, unplug it and have it checked by qualified personnel before further operation.
- Allow adequate air circulation to prevent internal heat build-up. Do not place the console on surfaces or near materials that may prevent proper ventilation.

- Do not install the console in a location near heat sources such as radiators or air ducts, and do not place console in direct sunlight. Isolate console from excessive dust, mechanical vibration or shock.
- Do not place the console on any unstable surfaces which may not properly support it. The console may fall resulting in injury and/or product damage. Use only with a cart, stand or table recommended by the manufacturer or sold with the console. A console and cart combination should be moved with care. Quick stops, excessive force and uneven surfaces may cause the console and cart combination to overturn.
- Grounding reliability can only be achieved when this equipment is connected to an equivalent receptacle marked "Hospital Grade."
- This product is for use only by qualified medical personnel trained for its use.
- **DANGER:** Risk of explosion if used in the presence of flammable anesthetics.
- If the unit has been subjected to a sudden temperature change, moisture may form on the metal inside the unit. If such a temperature change has occurred, allow the unit to achieve room temperature prior turning the unit on.
- Do not damage or modify the power cord. Damage to the power cord may cause a fire or shock hazard. When unplugging the power cord, always hold by the plug and remove it carefully.
- When not in use, unplug the console from the electrical outlet.
- When the console is not being used, store it in a location where it will not be damaged.
- Only ship the console in packaging that has been approved by NICO. Alternate packaging will not guarantee the protection of the device during shipping.
- Consult your distributor or NICO Customer Service if you can't correct a problem using this operator's manual.
- Precaution should be taken for those patients where increased risk or complications may be associated with the removal or biopsy of tissue based on the physician's judgment. Patients receiving anticoagulant therapy or who may have bleeding disorders may be at increased risk.
- The Myriad is not designed to coagulate tissue. Prior to resection of highly vascularized abnormalities, devascularization should be performed using alternate instruments.
- When operating the system in aspiration only mode, take care to avoid inadvertant contact with the side-kick foot pedal (to avoid unintentional enabling of cutting).
- The Myriad handpiece is designed to resect soft tissue. As such, care should be taken by the end user to limit resection to targeted tissue. Unintended resection of tissue may result in injury including but not limited to bleeding. Additionally, care should be taken during resection of vascular abnormalities as bleeding may occur.
- In accordance with the European REACH regulation and other environmental regulatory requirements, the Myriad console contains diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide, polyethylene glycol P-1,1,3,3-tetramethylbutylphenyl ether, Bisphenol-A and Methylene Chloride.

3 Myriad Console Introduction

3.1 Components (Included with Console)

- Foot Pedal
- Power Cord
- Nitrogen/Instrument Air Line
- Aspiration Line
- Cart


Refer to the Components and Accessories section in this Manual for a list of Myriad compatible components, including consumables. For additional information, please contact your local sales representative.

3.2 Console Functions and Symbols

Front Panel

NICO Myriad SPECTRA™ Console


NICO Myriad® and Myriad NOVUS® Consoles

1. System Power  - Toggle to turn power to the system on or off. I=ON, O=OFF
2. Foot Pedal  - Attachment site for the foot pedal cord to the console.
3. Prime Button  - This button enables the user to prime the handpiece, which includes simultaneous activation of cutting and aspiration (NOTE: maximum aspiration occurs during priming, regardless of the position of the aspiration control dial). Priming will occur as long as the prime button is depressed. The indicator light above the prime button illuminates while the handpiece is being primed. The cutter enable light, cutter active light and aspiration light will also illuminate while the handpiece is being primed.
4. Aspiration  - Attachment site for the aspiration line. The other end attaches to the top of the single-use collection canister. The aspiration line should always remain attached to the console, except when being replaced or the console is being shipped for service.

5. Aspiration Level Control Dial  - The user can increase or decrease the maximum aspiration level achieved at full depression of the foot pedal using this dial.
6. Aspiration Level Indicator - Digital display indicates aspiration level in non-descript units ranging from 0 to 25.
7. Handpiece  - Attachment site for the handpiece cord to the console. Indicator light above attachment site illuminates during cutting.
8. Cutter Enable Button  - This button enables/disables the cutting mode on the handpiece. The indicator light above this button indicates when cutting is enabled. An audible alert also indicates whether or not cutting is enabled. Two sequential tones indicate that cutting is enabled, while a single tone indicates that cutting is disabled.

Rear Panel

NICO Myriad® and Myriad NOVUS® Consoles

1. Power Cord Receptacle - Attachment site for detachable hospital grade power cord (IEC320C13).
2. Circuit Breaker - See “Circuit Breaker Reset” in Service and Maintenance section of this manual.
3. Rating Chart - Listing of the electrical specifications for the unit.

4. Nitrogen/Instrument Air Supply Input – Attachment site for the line. The other end of the line will attach to the source. The back of the console is fitted with the male DISS connector. The line is fitted with a female DISS connector at each end. If the OR supply has a different connection, adapters are available for purchase.

Symbols

	Handpiece - Attachment site for the handpiece cord to the console. Indicator light above attachment site illuminates during cutting.
	Cutter Enable Button - This button enables/disables the cutting mode on the handpiece. The indicator light above this button indicates when cutting is enabled. An audible alert also indicates whether or not cutting is enabled. Two sequential tones indicate that cutting is enabled, while a single tone indicates that cutting is disabled.
	System Power - Toggle to turn power to the system on or off. I=ON, O=OFF
	Foot Pedal - Attachment site for the foot pedal cord to the console.
	Prime Button - This button enables the user to prime the handpiece, which includes simultaneous activation of cutting and aspiration (NOTE: maximum aspiration occurs during priming, regardless of the position of the aspiration control dial). Priming will occur as long as the prime button is depressed. The indicator light above the prime button illuminates while the handpiece is being primed. The cutter enable light, cutter active light and aspiration light will also illuminate while the handpiece is being primed.
	Aspiration - Attachment site for the aspiration hose to the console. The other end of the hose attaches to the collection canister. The indicator light above attachment site illuminates during aspiration.
	Aspiration Level Control Dial - The user can increase or decrease the maximum aspiration level achieved at full depression of the foot pedal using this dial.
	Prescription Use Only
	Caution

	"Must Refer to Instructions for Use" Symbol
	Interference may occur in the vicinity of known RF transmitting devices and equipment marked with the following symbol.
	Nitrogen/Instrument Supply Input - Attachment site for the nitrogen or instrument air line to the console.
	Symbol for the Waste Electrical and Electronic Equipment Directive (WEEE Directive), European Community directive 2002/96/EC on waste electrical and electronic equipment (WEEE)
	Symbol indicating the catalog number for the device.
	Symbol indicating the serial number for the device.
	Symbol indicating the unique device identifier (UDI) for the device.
	Symbol indicating the location of an applied part. For the Myriad System, the applied part is the metal cannula on the handpiece.
	Symbol which specifies contact information for the European Authorized Representative.
	Symbol indicating the manufacturer of the equipment.
	Symbol indicating the date of manufacture of the equipment.
	Frequency Range
	AC Voltage
	Current Max
	Circuit Breaker
	Indicates item is a Medical Device

3.3 Console Connections

This section provides information for making all connections necessary to set up the console. The following connections will be covered:

- Power Cord
- Nitrogen/Instrument Air Line
- Aspiration Line
- Canister
- Foot Pedal
- Handpiece

Power Cord Connection

The power cord connects the console to the building power supply. Insert the hospital grade power cord into the power cord receptacle on the rear panel (Figure 1 & 2). The power cord is to be used for mains disconnection.

NOTE: Ensure that the power cord is in good condition. A damaged power cord poses an electrical shock hazard. When unplugging the unit, always grasp the plug and pull gently. NEVER pull on the cord.

Figure 1: Power Cord Connection

Figure 2: Power Cord Connection

Nitrogen / Instrument Air Line Connection

Supplied with the console are two approximately 18 ft (~5.5 m) lines with female DISS fittings on each end. One has a label and color code for Instrument Air, while the other is for Nitrogen. Choose proper color line for the hospital supply available. Insert one end of the line into the male fitting on the rear panel of the console (Figure 3 & 4). Tighten the connection by rotating the connection nut clockwise. The line should be hand tightened. Do not use tools to tighten the fittings. The opposite end of the hose should be connected to the hospital supply. Before attaching the line to the hospital supply, confirm the supply pressure does not exceed the recommended input supply pressure for the console. A DISS to Schrader adapter is available for Schrader connections. If the hospital supply connection is neither DISS nor Schrader, contact the biomedical engineering department, your distributor, or NICO Customer Service.

Figure 3: Nitrogen/Instrument Air Line Connection

Figure 4: Nitrogen/Instrument Air Line Connection

Aspiration Line Connection

The console will generate its own aspiration when activated by the foot pedal. During the initial console setup, a reusable flexible aspiration line is connected to the front of the console (Figure 5 & 6). To connect, firmly press the end of the tube into the fitting. The tube should insert approximately 0.5 inches (13 mm). The aspiration line should always remain attached to the console, except when being replaced or the console is being shipped for service.

Figure 5: Aspiration Line Connection

Figure 6: Aspiration Line Connection

The opposite end of the aspiration line includes a flexible suction connector that is attached to the vacuum port on top of the canister (Figure 7 & 8). Ensure that the lid and port caps are secure and sealed to avoid vacuum system leaks. Any leaks in the aspiration system will affect the performance of the Myriad.

Figure 7: Aspiration Line connected to the Canister Vacuum Port

Figure 8: Aspiration Line connected to the Canister Vacuum Port

Foot Pedal Connection

To connect the foot pedal to the console, insert the gray foot pedal connector into the gray receptacle on the console until it clicks into position (Figure 9 & 10) (NOTE: the foot pedal connector is keyed and will only connect to the console in one orientation). To remove the connector, pull back on outer sleeve with one hand while securing connector with other hand, then withdraw from console receptacle.

Figure 9: Foot Pedal Connection

Figure 10: Foot Pedal Connection

Handpiece Connection

To connect the handpiece to the console, insert the handpiece connector into the receptacle on the console until it clicks into position (Figure 11 & 12) (NOTE: the handpiece connector is keyed and will only connect to the console in one orientation). To remove the connector, pull back on the outer sleeve with one hand while securing connector in other hand, then withdraw from console and withdraw from console receptacle.

Figure 11: Handpiece Connection to Console


Figure 12: Handpiece Connection to Console

Figure 13 & 14 illustrate the handpiece tubing connected directly to the canister lid. This connection may be different if the end user chooses to collect tissue removed during the procedure. Additionally, components from the Automated Preservation System (APS) must be used to collect tissue. Additional information on the various APS setups can be found within the APS IFU. Section 6.7 contains a list of the components comprising the APS.

Figure 13: Handpiece Connection to Canister

Figure 14: Handpiece Connection to Canister

4 Using the Console

This section provides information for using the console.

4.1 System Power Up

After making all necessary connections, turn the system power on by depressing the switch to the "I" position (Figure 15 & 16). At this point, the indicator light next to the power switch will illuminate, other indicator lights on the user face will blink, and a single beep will be heard from the console. The aspiration indicator will read zero.

Figure 15: Console Power Up

Figure 16: Console Power Up

4.2 Priming the System

Before using the system, it is required that each handpiece be primed by aspirating fluid into the aperture, through the handpiece and tubing, and into the canister. This can be accomplished by depressing the prime button (Figure 17 & 18) while aspirating sterile fluid from a small bowl or from the reservoir provided in the tray in which the handpiece was packaged. This will ensure that the inside of the handpiece and tubing are lubricated and will reduce the likelihood of tissue occlusions.

Figure 17: System Prime

Figure 18: System Prime

4.3 Aspiration and Cutting

The console offers two modes of operation: (1) Aspiration only, and (2) Aspiration with Cutting.

This system is supplied with a two-stage foot pedal. The foot pedal consists of a large “main” foot pedal and a smaller “side-kick” foot pedal (Figure 19). When depressed, the large foot pedal will initiate either aspiration only or aspiration with cutting (depending on whether or not cutting has been enabled). The side-kick foot pedal enables the user to toggle between these modes of operation by either enabling or disabling cutting. The user may also toggle between these two modes by depressing the cutter enable button () on the console.

Figure 19: Foot Pedal

When cutting is enabled, depressing the large foot pedal will result in cutting action as well as aspiration. When cutting is disabled, depressing the large foot pedal will result in aspiration only. The level of aspiration delivered is proportional to the depth to which the foot pedal is depressed. Additionally, the aspiration dial on the console may be adjusted to control the maximum amount of aspiration delivered.

Refer to the Instructions for Use (IFU) provided with the handpiece for additional information concerning the handpiece.

To safely terminate operation of this equipment, cease depression of the foot pedal and power off the console.

5 Troubleshooting

This section provides guidance if problems are encountered when operating the console. Please consult the following information before contacting NICO to ensure the problem is not a misunderstanding of the operation of the system. If after reading this section the problem still cannot be resolved, please contact NICO for technical support.

Problem	Possible Cause	Action
No system response	Power cord not plugged in	Check cord. Make sure power cord is plugged in properly.
	Circuit breaker tripped	Reset circuit breaker on back of the console.
No aspiration/suction	Gas source not on	Check to make sure wall source or tank is on and tank is full
Handpiece body becomes warm/hot	Excessive tissue and/or blood in handpiece causing friction	Prime/flush handpiece vigorously with saline/ringers
Handpiece cutting slows down	Excessive tissue and/or blood in handpiece causing friction	Prime/flush handpiece vigorously with saline/ringers
Monitoring equipment in OR is affected by Myriad	Monitoring equipment not compliant with internationally recognized standards for electrical interference	Identify manufacturer and name of monitoring equipment, consult NICO Quality Department
Suction achieved via the foot pedal is less than when pressing the prime button.	Aspiration dial is set too low	Turn the aspiration dial up on the console
Handpiece not aspirating/cutting	Handpiece not plugged in	Check handpiece connection. Make sure handpiece is plugged into console.
	Aspiration line from console not connected to canister	Check console aspiration line. Make sure it is plugged into the canister.
	Aspiration line from handpiece not connected properly	Check the handpiece aspiration line. Make sure it is plugged into the canister or appropriate TPS accessory.
	Canister leak (crack in canister, crack in lid, uncapped port on lid, loose fit between canister and lid)	Check canister and lid for cracks, make sure lid is on tight, make sure all ports are capped
	Foot pedal not plugged in	Check foot pedal. Make sure foot pedal is plugged into console.
	Aspiration tubing clogged	Prime aspiration tubing. Make sure tubing is free from occlusions.
	Incorrect canister is used – volume is too large for effective vacuum delivery	NICO-specified canister MUST be used with the Myriad console
	Insufficient supply of gas to console – low tank	If a tank is being used, confirm that tank has adequate level. Readout on tank gauge must be at least 80 psi when foot pedal or prime button is fully depressed.
	Insufficient supply of gas to console – input pressure	Gas input pressure from wall or tank must be set at 100 psi.
	Dial setting on console is too low	Turn up aspiration dial on console
	Unvented scope during intraventricular case – all channels of scope are occupied	Explain that one channel of scope must remain open for venting so the handpiece will function properly.

5.1 Service and Maintenance

The NICO Myriad System is covered by a one-year warranty for defects, materials and workmanship. Additional extended service contracts are available. Contact your distributor or NICO Corporation directly for details.

NOTE: Repairs and adjustments inside the console are only to be performed by NICO. Customer attempts to repair/adjust the console will void the initial one-year warranty and/or any extended service contracts.

Circuit Breaker Reset

The circuit breaker is located on the rear panel of the console below the power cord receptacle.

1. Unplug console from the power source.
2. Disconnect power cord from receptacle on rear panel.
3. Allow 60 seconds for system to reset prior to proceeding to step 4.
4. Press the button on the circuit breaker to reset.
5. Plug the power cord to the console and into the power source.

Aspiration Line Replacement

To remove the aspiration line for replacement or shipping:

1. With thumb and forefinger, press the ring at the junction of the aspiration line and aspiration line fitting toward the NICO Myriad console.
2. While the ring is depressed, pull the aspiration line away from the console.

To replace the aspiration line:

1. Firmly press the end of the tube into the fitting. The tube should insert approximately 0.5 inches (13 mm).

Minimum Requirements for Hardware, IT Networks and Security Measures

There are no special requirements for hardware, IT networks or cyber-security. The device does not connect to a network and the firmware is inaccessible to the end-user.

5.2 Inquiries and Service

Direct inquiries to your distributor or to NICO Corporation:

NICO Corporation
Customer Service Department
250 East 96th Street, Suite 125
Indianapolis, IN 46240

Phone: +1-317-660-7118 ext.100

If your system requires service, contact your distributor or NICO Corporation Customer Service at +1-317-660-7118 ext.100. Your distributor or NICO Customer Service will provide you with all necessary information for returning and repairing your console.

5.3 Maintenance Schedule

Activity	Occurrence	Action
Wipe off console	As needed	Wipe off console with a clean towel moistened with a hospital approved solution
Wipe off foot pedal	As needed	Wipe off foot pedal with a clean towel moistened with a hospital approved solution
Inspect foot pedal cord	Annually	Check for cuts and damage to outside cover and strain reliefs
Inspect power cord	Annually	Check for cuts and damage to outside cover
Inspect Nitrogen/ Instrument Air line	Annually	Check for cuts and damage to outside cover
Test leakage current	Annually	Confirm that leakage current is acceptable per the requirements of the applicable hospital/facility
Inspect markings for legibility	Annually	View markings and confirm they are still legible (free from excessive fading, etc.)

NOTE: Of the activities listed above, only leakage current testing must be performed by hospital service personnel.

6 Specifications

6.1 NICO Myriad® and Myriad NOVUS® Consoles

Width.....	12.5 in. (318 mm)
Height.....	4.10 in. (104 mm)
Length.....	14.75 in. (375 mm)
Weight	16.50 Lbs. (7.0 kg)
Maximum Power	40 Watts
AC Voltage Range	100-240 V
Frequency Range.....	50-60 Hertz
Max Current	0.2 Amps
Circuit Breaker	2.0 Amps, 250 V
Nitrogen/Instrument Air Input Pressure Range.....	90-100 PSI (620-690 kPa)
Power Cord Length	15 feet (4.5 m) min.
Nitrogen/Instrument Air Line Length	18 feet (5.5 m) min.
Handpiece Cable Length.....	17 feet (5.2 m) min.

6.2 NICO Myriad SPECTRA™ Console

Width.....	12.3 in. (312 mm)
Height.....	4 in. (102 mm)
Length.....	14.4 in. (366 mm)
Weight	16.1 Lbs. (7.3 kg)
(Other Specifications are same as the Myriad and Myriad NOVUS Consoles)	

6.3 NICO Myriad® Foot Pedal

Width.....	4.75 in. (121 mm)
Height.....	3.75 in. (95 mm)
Length.....	11.875 in. (302 mm)
Weight.....	6.7 lbs. (3.0 kg)
Cord Length	17 feet (5.2 m)

6.4 Classification

Class I, Type BF Applied Part
Continuous Mode of Operation
Console: IP00
Foot Pedal: IP66

6.5 Console Equipment Operating Environmental Conditions

Temperature range within 15°C to 40°C
Relative humidity range within 30% to 75%
Atmospheric pressure range within 700 to 1060hPa
Operating altitude less than or equal to 3,000 meters above sea level

Storage Temperature : -29°C to 60°C
Storage Humidity : 30% to 85%
Transport Temperature : -29°C to 60°C
Transport Humidity : 30% to 85%

6.6 Disposal

Comply with all local codes when disposing of packaging, equipment or any other portion of this product.

6.7 Components & Accessories List

Name	Description	System Configuration
Capital Components		
Myriad Capital Kit	Cart Console Foot Pedal Nitrogen Line or Instrument Air Line	The console should be placed on the cart or suitable surface for use with the related device/components. When the cart is used, the foot pedal can be stored in the bracket the base of the cart.
Myriad NOVUS Capital Kit	Cart Console Myriad-LX Xenon Light Source Foot Pedal Nitrogen Line or Instrument Air Line	The console and light source should be placed on the cart or suitable surface for use with the related device/components.
Myriad SPECTRA Capital Kit	Cart Console Myriad SPECTRA Light Source Foot Pedals (Console and Light Source) Nitrogen Line or Instrument Air Line	When the cart is used, the foot pedal(s) can be stored in the bracket(s) at the base of the cart.
Name	Description	
Replacement Capital Components		
Nitrogen Line	18 ft. (5.5 m)	
Instrument Air Line	18 ft. (5.5 m)	
Power Cord (UK)	15 ft. (5.5 m)	
Aspiration Line	20 in. (508 mm)	
Required Consumable Components		
Canister	800 cc (mL) – NOTE: Use ONLY the canister supplied by NICO.	
Myriad Handpieces	Various Sizes	
Consumable Pack Options		
BrainPath Approach Packs	Various Sizes - includes BrainPath Sheath, Shepherd's Hook, Pre-Bent Myriad Handpiece, meniGLIDE, and SecureHold Retractors	
Illumination Packs	Various Sizes - Myriad Illumination Fiber and Myriad Handpiece Sleeve	
BrainPath Approach + Illumination Pack	Various - Combination of BrainPath Approach and Illumination Packs	
Automated Preservation System Packs	Various – Specimen Collector, Specimen Infusion Valve (SIV), Specimen Preserver, and Specimen Collector Identification Tags	
À La Carte and Optional Accessories		
Specimen Collector & Filter Element - Clamshell	Automated Preservation System (APS) accessories	
Replacement Filter Element - Clamshell		
Specimen Collector & Filter Element – Scoop		
Replacement Filter Element - Scoop		
Specimen Infusion Valve (SIV) with 0.5 mm Metering Line		
Specimen Infusion Valve (SIV) with 0.760 mm Metering Line		
Specimen Preserver		
Handpiece Universal Bending Tool	Handpiece Accessory	
Replacement PaediScope Handpiece Adapter	Handpiece Accessory	
Replacement LOTTA Handpiece Adapter	Handpiece Accessory	
Nitrogen Splitter	Non-Device Accessory – 3-way pneumatic	
DISS/Schrader Adapter	Non-Device Accessory – Non-Swivel Schrader x DISS Male	
DISS/Schrader Swivel Adapter	Non-Device Accessory – Swivel Schrader x DISS Male	
Schrader Quick Connect Adapter	Non- Device Accessory – Male DISS hand-tight / Female Schrader Adapter	

*Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution company, or customer service for the latest product offerings.

6.8 System Electromagnetic Emissions and Immunity Declarations

IEC 60601-1-2:2020 Table 1 Requirements			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the Equipment should assure that it is used in such an environment.			
Test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	Class A	The equipment is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the equipment or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A (Harmonics)	Class A (Harmonics)	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	10 Min	10 Min	

IEC 60601-1-2:2020 Table 4 Requirements			
Phenomenon Basic EMC standard or test method	Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radiated RF electromagnetic field IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications Equipment IEC 61000-4-3	See Table 9	See Table 9	
Power frequency (50 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity magnetic Fields IEC 61000-4-39	65 A/m at 134,2 kHz 7,5 A/m 13,56 MHz	65 A/m at 134,2 kHz 7,5 A/m 13,56 MHz	

IEC 60601-1-2:2020 Table 5 Requirements

Phenomenon Basic EMC standard or test method	Test level	Compliance level	Electromagnetic environment – guidance
Electrical fast transient/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency SIP/SOPS (if applicable) ± 1 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency SIP/SOPS (if applicable) ± 1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000- 4-5	$\pm 0,5$ kV, ± 1 kV line(s) to line $\pm 0,5$ kV, ± 1 kV, ± 2 kV line(s) to ground	$\pm 0,5$ kV, ± 1 kV line(s) to line $\pm 0,5$ kV, ± 1 kV, ± 2 kV line(s) to ground	Mains power quality should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 UT = 230 Vac	Voltage Dips: 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle And 70 % UT; 25/30 cycles Single phase: at 0° Voltage Interruptions: 0 % UT; 250/300 cycle	Voltage Dips: 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle And 70 % UT; 25/30 cycles Single phase: at 0° Voltage Interruptions: 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment requires continued operation during power mains interruptions, it is recommended that the equipment be powered from an uninterruptible power supply or a battery.

IEC 60601-1-2:2020 Table 9 Requirements

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation	Immunity Test Level
385	380 to 390	TRTRA 400	Pulse Modulation ^{b)} 18 Hz	27
450	430 to 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	28
710	704 to 787	LTE Band 13, 17	Pulse Modulation ^{b)} 217 Hz	9
745				
780				
810	800 to 960	GMS 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation ^{b)} 18 Hz	28
870				
930				
1720	1700 to 1990	GMS 1800; CDMA 1900; GMS 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation ^{b)} 217 Hz	28
1845				
1970				
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation ^{b)} 217 Hz	28
5240	5100 to 5800	WLAN 802.11 a/n	Pulse Modulation ^{b)} 217 Hz	9
5500				
5785				

If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.

7 General Information

Special Storage or Handling Conditions

There are no special storage or handling conditions.

Device Related Serious Incident Reporting

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at +1-317-660-7118 ext.100.

Major Risks from Risk Management Activities

- Electric Shock
- Fire/Explosion
- Damage to Healthy Neurological Tissue
- Excessive Bleeding
- Infection
- Transmission of Disease
- Loss of Function of Equipment in Operating Room due to EMI Emissions
- Tissue Irritation
- Pyrogenic Response (Fever)

The published literature and NICO's post-market complaint processes continue to support very low occurrence or absence of these risks. All residual risks associated with the Myriad continue to be acceptable when weighed against the potential clinical benefits.

Clinical Benefits:

- Targeted tissue resection
- Favorable safety profile

Summary of Safety & Clinical Performance

Where applicable, a summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

Basic-UDIs:

Capital Equipment	081674402MY3FY
Handpieces	081674402MY1FU
Illumination Accessories	081674402MY2FW
Preservation & Other Accessories	081674402MY4G2
	081674402MY5G4

For More Information

For more information on the NICO Myriad, contact your NICO representative or NICO Customer Service at +1-317-660-7118 ext.100. Additional information may be found at www.niconeuro.com. Additional information may be found at www.niconeuro.com, including electronic versions of the labelling: www.niconeuro.com/labeling

Patent: www.niconeuro.com/patents

© 2025 NICO Corporation. All Rights Reserved. Printed in the USA.

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© 2025 NICO Corporation. Reservados todos los derechos. Impreso en los Estados Unidos.