

meniGLIDE® Product Information



Provides access through the meninges for intracranial and spine procedures.

The meniGLIDE is designed for **control** and **safety** to open the meninges while protecting collateral tissue in preparation for a surgical procedure.


Product Features

1

Malleable, low-profile shaft for enhanced visualization

2

Atraumatic dissecting tip and footplate designed to protect underlying structures



*This image has been magnified for more clear visualization of the product features.


.5 cm length
dissecting blade

1cm length
footplate

3.4 inch ergonomic handle

Product Information

- Product Number: NN-8040
- Sterile / Single-use
- May be purchased as an ancillary product (5 per box)



Click or scan to view the meniGLIDE achieving access



NICO and NICO meniGLIDE are registered trademarks of NICO Corporation in the United States and other countries. All rights reserved. Patents pending worldwide. Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution company, or customer service for the latest product offerings.

